
# 20-40X ZOOM MICROSCOPE

- 1) Eyepiece
- ② Zoom Wheel
- (3) LED On/Off Button
- 4 Focusing Lever
- (5) Battery Insulate
- (6) Specimen Stage
- (7) Battery Door
- (8) Battery Door Screw
- (9) Specimen Slides Covers
- 10 Specimen Slide



## **HOW TO USE**

#### Getting Ready


Remove the Battery Insulate (5) by pulling it outward.


#### Observing with Specimen Stage (6)

- Place the object of observation onto the Specimen Slide (10) and place a Specimen Slide Cover (9) on top of the object. Insert the Specimen Slide into the Specimen Stage.
- Place the Microscope into the Specimen
  Stage and press the LED On/Off Button
  (3) and observe through the Eyepiece (1).
- 3. Turn the Zoom Wheel to adjust the desired degree of magnification (from 20X to 40X).
- 4. Slide the Focusing Lever (4) gradually until the image is clear and sharp.


#### Observing without Specimen Stage

- Place the microscope directly on top of the object of observation and press the LED On/Off Button and observe through the Eyepiece.
- Turn the Zoom Wheel to adjust the desired degree of magnification (from 20X to 40X).
- 3. Slide the Focusing Lever (4) gradually until the image is clear and sharp.





### **BATTERY INSTALLATION**



If the LED light becomes dim, it means that it's time to replace the batteries.

- Unscrew the Battery Door Screw (8) with a screwdriver and remove the Battery Door (7) in the direction shown by the arrow.
- 2. Remove the used-up batteries and insert 3 pcs of fresh LR44 button cells into the battery compartment in the right polarity.
- 3. Close the Battery Door and tighten the Battery Door Screw with a screwdriver.