**OPTIKA**

M I C R O S C O P E S  
I T A L Y

Ver. 5.0.0

# **B-20 & B-20Solar**

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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## INDEX

<b>1.0 DESCRIPTION</b>	<b>page 4</b>
<b>2.0 INTRODUCTION</b>	<b>page 6</b>
<b>3.0 UNPACKING AND ASSEMBLY</b>	<b>page 6</b>
<b>4.0 USING THE MICROSCOPE</b>	<b>page 6</b>
<b>5.0 MAINTENANCE</b>	<b>page 7</b>
<b>6.0 TECHNICAL SPECIFICATIONS</b>	<b>page 8</b>
<b>7.0 RECYCLING AND RECOVERY</b>	<b>page 8</b>

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

## Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

## Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

## Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

### DANGER

The instructions on this symbol to avoid possible severe personal injuries.

### WARNING

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

### WARNING

Possibility of electric shock.

### HOT!

Attention: high temperature surfaces. Avoid direct contact.


### NOTE

Technical notes or usage tips.

## 1.0 DESCRIPTION

**B-20: FRONT VIEW**

**B-20: REAR VIEW**

## 2.0 INTRODUCION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

## 3.0 UNPACKING AND ASSEMBLY

- 3.1 The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Lift the microscope from the container and put it on a stable desk.
- 3.2 Insert the eyepiece into the eye tube.

**Connect the mains plug into the socket at the base**

### WARNING

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

### WARNING

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

### WARNING

#### AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

## 4.0 USING THE MICROSCOPE

### 4.1 Place the specimen on the stage

Turn the observation head to a comfortable position for observation. Lock the specimen slide on the mechanical stage using the slide clip. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

### 4.2 Lamp settings

The microscope is fitted with white LED illuminator. Before turning on the light, read the section 5.3 about electrical safety precautions. Insert the plug of the cable into the power socket and turn on the

switch on the back of the main body. Turn the brightness adjustment knob to a brightness suitable for observation.

### 4.3 Focus adjustment

Rotate the coarse focusing knob to bring the slide into focus using an objective with a low magnification. Adjust the fine focusing knob to get it in focus, turn the nosepiece to choose the objective you need.

## 5.0 MAINTENANCE

### 5.1 Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### 5.2 To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### 5.3 Electrical safety precautions

- Before plugging in the power cord with the supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Do not turn the power on and off, off and on immediately as this will shorten the life span of the bulb and may cause damage to the electrical system.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### 5.4 Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.

Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.

- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

### 5.5 If you need to send the microscope to Optika for maintenance, please use the original packaging


## 6.0 TECHNICAL SPECIFICATIONS

<b>HEAD:</b>	Monocular, 360° rotating, 45° inclined
<b>EYEPIECES:</b>	WF10X/16mm
<b>NOSEPIECE:</b>	Triple, reversed
<b>OBJECTIVES:</b>	Achromatic 4x, 10x, 40x
<b>WORKING STAGE:</b>	Rotating round stage, diameter: 90 mm, moving range 5 mm, slide clips
<b>FOCUSING SYSTEM:</b>	Coarse and fine with different axis, with focus tension adjustment
<b>ILLUMINATION:</b>	White LED, non-rechargeable, with brightness control
<b>MAINS POWER:</b>	220V / 50Hz B-20SOLAR: solar cell
<b>FUSE:</b>	2A

## 7.0 RECOVERY AND RECYCLING

Art.13 DLsg 25 july 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**New York Microscope Company Inc**