**B-380 Series**

# INSTRUCTION MANUAL

<b>Model</b>
B-383MET

Version: 1  
Issued: 18, 05, 2014

---

## Table of Contents

**Warning**

**Symbols and conventions**

**Safety Information**

**Intended use**

**Package Contents**

**Unpacking**

*Product Specifications*

*Overview*

**Using the microscope**

**Maintenance**

**Troubleshooting**

**User replaceable accessories and spare parts**

**Equipment disposal**

## Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance.

It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

## Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

## Package Contents

DESCRIPTION	QUANTITY
Microscope stand with nosepiece, stage, condenser	1
Optical head (trinocular)	1
Objective PLAN IOS MET 4x	1
Objective PLAN IOS MET 10x	1
Objective PLAN IOS MET 20x	1
Objective PLAN IOS MET 50X	1
Eyepiece WF10X/20mm	2
Metallographic attachment	1
Polarizer & Analyzer filters	1
Glass insert for sample support	1
Dust cover	1
Power supply output 6Vdc	1

---

## Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.


Place the metallographic attachment onto the top of the stand and tighten the lock-screw. Then place the observation head onto the top of the attachment and tighten the other lock-screw. Insert the eyepieces into the eye tubes and fix them with the tiny screws on the side of the eye tubes. Remove the plastic protection film from the microscope stage.


Connect the cable from the epillumination LED to the upper jack on the rear, and the provided +6Vdc power supply to the lower jack (marked with "6V" label).

## Product specifications

<b>Head:</b>	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
<b>Eyepieces:</b>	WF10X/20mm.
<b>Nosepiece:</b>	5-position reversed revolving nosepiece. Ball bearing rotation.
<b>Objectives:</b>	PLAN IOS MET 4x/0.10, 10x/0.25, 20x/0.40, 50x/0.75.
<b>Focusing system:</b>	Coaxial coarse and fine.
<b>Stage:</b>	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm. Belt-drive in X direction. Glass insert for sample support.
<b>Condenser:</b>	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
<b>Illumination:</b>	Transmitted light: X-LED <sup>3</sup> , with manual brightness control.
<b>Power supply:</b>	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A.
<b>Metallurgical Attachment:</b>	High-power white X-LED <sup>3</sup> ; field and aperture diaphragms.
<b>Filters:</b>	Blue, polarizer & analyzer.

## Overview

---

## Using the microscope

### Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

### Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage. A glass plate is provided in order to support metallurgical samples.

### LED lamp settings – Transmitted illumination

The microscope is fitted with high brightness white LED illuminator. Press the selector switch on the back of the main body to I position. Turn the brightness adjustment knob to a brightness suitable for observation.

### LED lamp settings – Epillumination

The microscope is fitted with high brightness white LED illuminator. Press the selector switch on the back of the main body to II position. Turn the brightness adjustment knob to a brightness suitable for observation.

### Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

### Focus and diopter adjustment

Turn the dioptic adjustment ring on the left eye tube to the “0” position. Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus using an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptic adjustment ring until a sharp and clear image can be observed also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

### Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

### Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution. There are two aperture diaphragms on the microscope: for transmitted light and for epillumination light.

### Adjust the field diaphragm for epillumination

Move the lever on the metallographic attachment to fully close the field diaphragm. Then slowly open it and stop when the whole field of view is illuminated. This condition prevents the light that falls outside the field of view from being scattered and producing noise.

### Polarisation and colour filters

The blue, yellow and frosted glass filters can be inserted in the slot in the metallographic attachment, just in front of the Led house. In the same slot you can insert the polarizer filter. The analyzer filter can be inserted in the slot below the optical head, after removing the protection cap.

Turn the graduated wheel on the analyzer filter to perform polarisation analysis of your sample.

## Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.  
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>LIGHT DOESN'T TURN ON</b>	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
<b>IMAGE CANNOT BE SEEN OR IS DARK.</b>	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
<b>IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.</b>	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination


## User replaceable accessories and spare parts

CAT. NO.	DESCRIPTION
M-160	Eyepiece WF10x/20mm.
M-161	Eyepiece WF15x.
M-162	Eyepiece WF20x.
M-163	Micrometer eyepiece WF10x/20mm.
M-005	26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm.
M-144	Objective 4x/0,10 E-PLAN IOS.
M-145	Objective 10x/0,25 E-PLAN IOS.
M-146	Objective 20x/0,40 E-PLAN IOS.
M-147	Objective 40x/0,65 E-PLAN IOS
M-149	Objective 60x/0,80 E-PLAN IOS.
M-148	Objective 100x/1,25 E-PLAN IOS (Oil).
M-337	Objective PLAN IOS MET 4x/0.10.
M-338	Objective PLAN IOS MET 10x/0.25.
M-339	Objective PLAN IOS MET 20x/0.40.
M-335	Objective PLAN IOS MET 50x/0.75.
M-173	Photo adapter for reflex camera.
M-114	CCD camera adapter.
M-069	Solar battery pack.

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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## **OPTIKA® S.r.l.**