OPERATION MANUAL

GUIDA UTENTE

MANUAL DE INSTRUCCIONES

MANUEL D'INSTRUCTIONS

1.0 DESCRIPTION	page 4
2.0 INTRODUCTION	page 6
3.0 UNPACKING AND ASSEMBLY	page 6
4.0 USING THE MICROSCOPE	page 8
5.0 MAINTENANCE	page 12
6.0 ELECTRICAL SPECIFICATIONS	page 13
7.0 RECYCLING AND RECOVERY	page 14

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

DANGER

The instructions on this symbol to avoid possible severe personal injuries.

WARNING

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

WARNING

Possibility of electric shock.

HOT!

Attention: high temperature surfaces. Avoid direct contact.


NOTE


Technical notes or usage tips.

1.0 DESCRIPTION

SPECIFIC PARTS FOR BRIGHTFIELD ILLUMINATION:

2.0 INTRODUCCION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily classroom and laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.


3.0 UNPACKING AND ASSEMBLY

- 3.1** The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.
- 3.2** The objectives are located in individual protective vials. Remove the objectives from the vials and insert them into the microscope nosepiece in the order from the lowest magnification to the highest, in a clockwise direction from the rear.
- 3.3** Epifluorescence parts (epiilluminator with filters block, power supply and mercury lamp housing) are housed in a hard case. Remove the protection caps from the top of the main B-600 body and from the epiilluminator and place the epiilluminator on the top of the microscope body. Tighten the screw on the right side using the provided allen wrench.

EPIILLUMINATOR
SCREW

- 3.4** Insert the mercury lamp house on the back of the epiilluminator and tighten the screw using the provided allen wrench.


LAMP SCREW

3.6 Insert the eyepieces into the eye tubes.

3.7 Fix the orange protection screen to the front of the epiilluminator with the two provided screws.

3.8 Mercury lamp insertion and replacement

Disconnect any cable from the epiilluminator. Loosen the lock-screw on the top of the lamp housing. The lamp holder will come out with the door.

Remove the old lamp or the plastic placeholder. Use protection glasses as, if handled carelessly, there is a risk that the lamp may shatter and eject glass splinters. Thus, without touching the glass part of the lamp, take out the fluorescence lamp from its box, taking care not to exert stress on any part of the lamp. Put one end inside the flexible metal holder (observe that it fits only in one direction) and then in the fixed one. Tighten the screws just enough to make it stay in place. If you tighten them too much the lamp may shatter. Put the lamp door with lamp assembly back in place and connect the cable on the side of the epiilluminator to the power supply. Connect the mains cable of the power supply. Turn on the power supply using the power switch. Wait until the current display shows about 5.0 A. If the current falls below 4.5A, replace the lamp. You should now wait for at least 10 minutes before aligning the lamp.

On the life time display you can read for how many hours the lamp has been in use. Replace the lamp every 500 hours, even if it seems to work correctly, as otherwise it may shatter. Take care to let the lamp cool down at least 30 minutes before opening the lamp housing. After having replaced the lamp, press the reset button to reset the life time count.

3.0 UNPACKING AND ASSEMBLY

Connect the mains plug into the socket at the base

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

4.0 USING THE MICROSCOPE

4.1 Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

4.2 Place the specimen on the stage


Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

4.3 Epifluorescence use: Aligning the Hg lamp

- Move the filter selector of the filters block on the desired position, corresponding to an excitation wavelength.
- Insert a specimen slide and focus with 4X objective.
- Close completely the Field Diaphragm by pulling down the lever marked "FS". You should see the diaphragm in the field of view. Act on the pair of knobs on both sides in order to center the diaphragm.

- Fully open the Field Diaphragm. Now switch to a greater magnification objective and center the lamp in order to maximize the light intensity (acting on the lamp centering knobs).
- Adjust the rotating lever for optimization of the lamp focusing.

- The lamp house allows also a fine tuning of the orientation of the mirror located at the back of the lamp. Use the provided allen wrench to change the orientation.

For more information about the lamp, see section 5.

- 4.4** Always keep the iris diaphragm completely open when working with fluorescence light.
- 4.5 Brightfield use: Lamp settings**
The microscope is fitted with a 3.5W X-LED³ illuminator. Insert the plug of the cable into the power socket and turn on the switch on the back of the main body. Turn the brightness adjustment knob to a brightness suitable for observation.
- 4.6 Adjust interpupillary distance**
Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. After finishing the interpupillary adjustment, turn the dioptic adjustment ring on the left eye tube to the “0” position at the graduation line on the eyetube.
- 4.7 Focus and diopter adjustment**
Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus with an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptic adjustment ring to a sharp and clear image also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.
- 4.8 Brightfield use: Condenser adjustment**
Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser-centring screws can be used to centre the condenser.
- 4.9 Brightfield use: Set the numerical aperture**
Adjust the aperture of the iris diaphragm located under the condenser, to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.
- 4.10 Brightfield use: Field diaphragm**
Adjust the field diaphragm to align the system and to control stray light.
- 4.11 Brightfield use: Additional filters**
The chromatic correction filter (blue) and the monochromator filter (green) can be inserted into the filter holder for definition of specimen parts. The use of the green filter is recommended for phase contrast observations.
- 4.12 Brightfield use: Using the polarizer (optional)**
The polarizing set consists of two filters; the analyser and the polarizer. The analyser is placed in the analyser holder (you have to remove the orange protection filter in order to access the holder) and the polarizer on top of the illuminator.
- 4.13 Video capturing (optional)**
B-500TiFL is provided with trinocular observation head. It can be connected to cameras via a photo/video adaptor, for photo and video capturing. Before taking a picture or filming video, if necessary, obscure with a dark cloth both the camera viewfinder and the eyepieces and pull out the light path selector lever. Please refer to the adaptor and camera manuals for further details.



4.14 Using the Phase Contrast set (optional)

The phase contrast set has 4 objectives for phase contrast and one condenser, with 4 rings of phase, and one position for brightfield. To use the phase contrast, remove the Abbe condenser of the microscope and replace it with the phase contrast condenser.

If, during use, the phase contrast diaphragm rings (anuli) are dislodged, they can be realigned rotating the alignment screws placed on two sides of each ring, using the enclosed tools. In particular:

-Select an objective corresponding to the ring to be realigned, and focus on a slide under brightfield illumination.

- The selector of the condenser has 5 selectable positions:
 - 0 → to use when it is necessary to be used the instrument in brightfield
 - 10 → to use with 10x phase contrast objective
 - 20 → to use with 20x phase contrast objective
 - 40 → to use with 40x phase contrast objective
 - 100 → to use with 100x phase contrast objective
- Insert the phase ring to be realigned.
- Remove an eyepiece, and replace it with the centring telescope.
- Focus the centring telescope, without touching the microscope focus, in order to focus the phase ring and the phase plate simultaneously. In particular, the phase plate is visible as a dark ring, and the phase ring as a bright ring divided into sectors.
- Turning the two phase ring centring screws, position the phase ring (bright ring) so that it is completely contained in the phase plate profile (dark ring). (See figure.)
- Remove the telescope, and replace the eyepiece.

4.15 Darkfield condenser (optional)

For darkfield observation, the Abbe condenser in standard equipment must be replaced with a special condenser for darkfield, showed in the picture.

5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing).
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.2 Do NOT!


- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

5.3 Safety precautions - fluorescence mercury lamp

- Always turn on the lamp with the current adjusted to maximum (5 A) and then wait for a minimum of 10 minutes before using the microscope.
- If the maximum current goes down below 4.5 A the lamp needs to be changed.
- After the lamp has been switched off it takes about 30 minutes before it is cool enough to be ready to be used again.
- Be sure the lamp is cool before you cover the microscope with the dust-cover after use.
- The lamp has a lifetime of 500 hours. In order to avoid explosive failure, replace the lamp after 500 hours, even if it seems to work correctly.
- The lamp light contains ultraviolet radiation that is harmful for the eyes and skin. Always look at the lamp light through the orange shield or through the dark alignment window on the filter block.

5.4 How to change the fluorescence mercury lamp

Disconnect any cable from the epiilluminator. Loosen the lock-screw on the top of the lamp housing. The lamp holder will come out with the door.

Remove the old lamp or the plastic placeholder. Use protection glasses as, if handled carelessly, there is a risk that the lamp may shatter and eject glass splinters. Thus, without touching the glass part of the lamp, take out the fluorescence lamp from its box, taking care not to exert stress on any part of the lamp. Put one end inside the flexible metal holder (observe that it fits only in one direction) and then in the fixed one. Tighten the screws just enough to make it stay in place. If you tighten them too much the lamp may shatter. Put the lamp door with lamp assembly back in place and connect the cable on the side of the epilluminator to the power supply. Connect the mains cable of the power supply. Turn on the power supply using the power switch. Wait until the current display shows about 5.0 A. If the current falls below 4.5A, replace the lamp. You should now wait for at least 10 minutes before aligning the lamp.

On the life time display you can read for how many hours the lamp has been in use. Replace the lamp every 500 hours, even if it seems to work correctly, as otherwise it may shatter. Take care to let the lamp cool down at least 30 minutes before opening the lamp housing. After having replaced the lamp, press the reset button to reset the life time count.

5.5 If you need to send the microscope to Optika for maintenance, please use the original packaging.

6.0 ELECTRICAL SPECIFICATIONS

Input Voltage:	100 - 240 Vac, 50/60 Hz
Brightfield Lamp:	3.5W X-LED ³
Fluorescence Supply:	100 - 240 Vac, 50/60 Hz
Fluorescence Lamp:	HBO 100 W high pressure mercury bulb

7.0 RECOVERY AND RECYCLING

Art.13 Dlsg 25 july 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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New York Microscope Company Inc