HDMI Cameras – H Series

INSTRUCTION MANUAL

Model
C-HP4

Ver. 1.0 2021

Table of Contents

1.	Warning	3
2.	Symbols and conventions	3
3.	Safety Information	3
4.	Intended use	3
5.	C-HP4 Basic features	4
5.1	Package content	4
5.2	Quick Instructions for C-HP4 camera	5
5.2.1	HDMI Mode	5
5.2.2	USB Mode	5
6.	C-HP4 User Interface (UI) and its Functions	6
6.1	The Camera Control Panel	7
6.2	The Measurement Toolbar	8
6.3	Camera calibration	9
6.4	Icons and Functions of the Camera Control Toolbar	11
6.4.1	Setting > Measurement	11
6.4.2	Setting > Magnification	12
6.4.3	Setting > Image Format	12
6.4.4	Setting > Video	13
6.4.5	Setting > Storage	14
6.4.6	Setting > Files	14
6.4.7	Setting > Language	15
6.3.8	Setting > Miscellaneous	16
7.	Micrometric Slide M-005	17
	Equipment disposal	18

1. Warning

This device is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. C-HP4 Basic features

C-HP4 is a multiple interfaces (HDMI + USB2.0 + SD card) CMOS camera.
HDMI and USB2.0 are used as the data transfer interface to HDMI display or computer.

- For HDMI output, the *Camera Control Panel + Measurement Toolbar* and *Camera Control Toolbar* are overlaid on the HDMI screen when the mouse move to the related region.
In this case, the USB mouse can be used to set the camera, browse and compare the captured image, play the video and perform the measurement.
- For USB output, connect computer with USB connection to transfer video in real time.
From here the ProView software can control each camera function.

5.1 Package content

- | | |
|------------------|---------------------|
| ① Camera | ⑤ Micrometric slide |
| ② HDMI cable | ⑥ SD card 16GB |
| ③ USB2.0 cable | ⑦ Power supply |
| ④ Wireless mouse | |

5.2 Quick Instructions for C-HP4 camera

Before turning on the camera, please connect the standard C-mount camera to the camera adapter and insert it into the microscope's trinocular port.

Figure 1 - C-HP4 Rear Panel

- | | |
|-----------------------|---------------------|
| ① USB mouse connector | ⑤ SD card slot |
| ② USB2.0 connector | ⑥ Power supply plug |
| ③ HDMI connector | ⑦ LED indicator |
| ④ ON-OFF button | |

You can use the C-HP4 camera in 2 different ways. Each application requires different hardware environment.

5.2.1 HDMI Mode

1. Plug the HDMI cable into the HDMI port ③ to connect the C-HP4 camera to HDMI display.
2. Plug the USB mouse into USB Mouse port ① to get control of the camera by using built-in software LiteView.
3. Plug 12V/1A power adapter into DC 12V slot ⑥ to supply power for the camera. The LED Indicator ⑦ will turn into red.
4. Insert SD card into SD Card slot ⑤ for saving captured images and recorded videos.
5. Press ON/OFF Button ④ to turn on the camera. The LED indicator will turn into blue.
6. Move mouse cursor to the left side of the video window, a *Camera Control Panel* will appear. It includes Manual/ Automatic Exposure, White Balance, Sharpness and other functions, please refer to 6.1 for details.
7. Move mouse cursor to the bottom of the video window and a *Camera Control Toolbar* will appear. Operations like Zoom In, Zoom Out, Flip, Freeze, Cross Line, Comparison and so on can be achieved. Please refer to 6.3 for details.
8. Move mouse cursor to the upper side of the video window, a *Measurement Toolbar* with calibration and other measurement tools will appear, please refer to 6.2 for details. The measurement data can be output with *.CSV format.

5.2.2 USB Mode

1. Plug the USB cable provided with the camera into USB2.0 ② to connect the C-HP4 camera to the computer.
2. Plug 12V/1A power adapter into DC 12V slot ⑥ to supply power for the camera. The LED Indicator ⑦ will turn into red.
3. Press ON/OFF button ④ to turn on the camera. The LED indicator will turn into blue.
4. Start the ProView software.
5. Clicking the camera name in the camera list starts the live video.

- **NOTE: When the USB cable and the mouse are plugged into the camera at the same time, the USB cable has priority and the mouse is not available; when the USB cable is unplugged, the mouse can be used normally.**

6. C-HP4 User Interface (UI) and its Functions

The C-HP4 UI shown in Fig. 2 includes a *Camera Control Panel* on the left side of the video window, a *Measurement Toolbar* on the upper side of the video window, a *Camera Control Toolbar* on the bottom of the video window.



Figure 2 - Camera Control User Interface

1. Move the mouse cursor to the left side of the video window: the *Camera Control Panel* will pop up automatically.
2. Move the mouse cursor to the upper side of the video window: a *Measurement Toolbar* will pop up for the calibration and measurement operations.
 - When user left-clicks the *Float/Fixed* button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if users move mouse cursor to left side of the video windows.
 - Only when user left-clicks the  button on the *Measurement Toolbar* to exit from measuring procedure it will be possible to perform other operations on the *Camera Control Panel*, or *Camera Control Toolbar*.
 - During the measuring process, when a specific measuring object is selected an *Object Location & Attributes Control Bar*  will appear for changing location and properties of the selected objects.
3. Move the mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically.

6.1 The Camera Control Panel

The *Camera Control Panel* controls the camera to achieve the best image quality according to the specific applications. It will pop up automatically when mouse cursor is moved to the left side of the video window (in measurement status, the Camera Control Panel will not pop up. Only when measurement process is terminated the Camera Control Panel will pop up by moving mouse cursor to the left side of the video window). Left-clicking button  to achieve *Display/ Auto Hide* switch of the Camera Control Panel.

Control Panel	Function	Description
	Snap	Capture image from the current video window
	Record	Record video from the current video window
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	Exposure Compensation	Available when Auto Exposure is checked. Slide to left or right to adjust Exposure Compensation according to the current video brightness to achieve proper brightness value
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	Auto White Balance	White Balance adjustment according to the window video every time the button is clicked
	Manual White Balance	Slide Red or Blue cursors to manually set the video White Balance
	ROI White Balance	When the ROI region is moved, white Balance can be adjusted according to content inside the ROI region.
	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video
	Green	Slide to left or right to decrease or increase the proportion of Green in RGB on video
	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video
	Sharpness	Adjust Sharpness level of the video window
	Denoise	Slide left or right to denoise the image
	Saturation	Adjust Saturation level of the video window
	Gamma	Adjust Gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma
	Contrast	Adjust Contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast
	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check AC(50HZ) to eliminate flickering "strap" caused by 50Hz illumination
	AC(60HZ)	Check AC(60HZ) to eliminate flickering "strap" caused by 60Hz illumination
Default	Restore all the settings in the Camera Control Panel to default values	

6.2 The Measurement Toolbar

The *Measurement Toolbar* will pop up when moving mouse cursor to any place near the upper side of the video window.

Figure 3 - The Measurement Toolbar

Icon	Function
	<i>Float/ Fix</i> switch of the <i>Measurement Toolbar</i>
<input checked="" type="checkbox"/> Visible	Define measuring object in <i>Show up/ Hide</i> mode
<input type="text" value="Nanometer (nm)"/>	Select the desired <i>Measurement Unit</i>
<input type="text" value="4X"/>	Select <i>Magnification for Measurement after Calibration</i>
	<i>Object Select</i>
	<i>Angle</i>
	<i>4 Points Angle</i>
	<i>Point</i>
	<i>Arbitrary Line</i>
	<i>3 Points Line</i>
	<i>Horizontal Line</i>
	<i>Vertical Line</i>
	<i>3 Points Vertical line</i>
	<i>Parallel</i>
	<i>Rectangle</i>
	<i>Ellipse</i>
	<i>5 Points Ellipse</i>
	<i>Circle</i>
	<i>3 Points Circle</i>
	<i>Annulus</i>
	<i>Two Circles and Center Distance</i>
	<i>3 Points Two Circles and Center Distance</i>
	<i>Arc</i>
	<i>Text</i>

	Polygon
	Curve
	Scale Bar
	Arrow
	Execute <i>Calibration</i> to determine the corresponding relation between magnification and resolution, which will establish the corresponding relationship between measurement unit and the sensor pixel size. <i>Calibration</i> needs to be done with the help of a micrometer. For detailed steps of carrying out <i>Calibration</i> please refer to chapter 6.3.
	Export the <i>Measurement</i> to CSV file(*.csv)
	<i>Measurement Setup</i>
	Delete All the Measurement Objects
	Exit from <i>Measurement</i> mode
	When the measurement ends, left-click on a single measuring object and the <i>Object Location & Properties Control Bar</i> will show up. User could move the object by dragging the object with the mouse. But more accurate movement could be done with the control bar. The icons on the control bar mean <i>Move Left, Move Right, Move Up, Move Down, Color Adjustment</i> and <i>Delete</i> .


Note:


1. When user left-clicks *Display/Hide* button  on the *Measurement Toolbar*, it will be fixed. In this case the *Camera Control Panel* will not pop up automatically even if moving the mouse cursor to the left edge of the video window. Only when user left-click the  button on the *Measurement Toolbar* to exit from the measurement mode it will be possible to perform other operations on the *Camera Control Panel* or on the *Camera Control Toolbar*.
2. When a specific item is selected during the measurement process, the *Object Location & Attributes Control Bar*  will appear for changing the object location and properties of the selected objects.

6.3 Camera calibration


The camera must be calibrated before any measurements can be taken. To do this it is necessary to use the micrometer slide supplied with the camera.

Detailed procedure is the following:


1. Insert the objective with the lowest magnification available and put the stage micrometer on the stage.
2. Click on the “*Calibrate*”  button in the Measurement toolbar.
3. The following dialog box appears:

4. In the *Magnification* field enter the current magnification
 - Some magnification values (default ones 4x, 10x, 20x, 40x and 100x) are preset in the drop-down list.
 - If user must insert a different magnification value, double click in the *Magnification* field. The following dialog box appears:

- From here user can enter any desired magnification value.
5. Focus the ruler under the microscope.
 - Along with the previous dialog box a red line appears on the image. This line is the calibration line.

6. Align the two ends of the red line with the ruler scale (try to use as many notches as possible to keep the calibration accuracy).
7. Set the desired unit (micrometer, millimeter or nanometer) in the *Calibration* dialog box.
8. Insert the length of the drawn line in the *Actual Length* field (using the same procedure described above).
 - Micrometric slide has two rulers (see chapter 11). One is 1 mm/100 (for biological microscopes) and the other is 10mm/100 (for stereomicroscopes).
 - When using the ruler 1mm/100 the interval between two long lines is 100 μm , while using the ruler 10mm/100 the interval between two long lines is 1000 μm .
9. If everything is fine, click *OK* to end the calibration. The new magnification (the number you entered in the *Magnification* edit box) will be available in the *Magnification* drop-down list.
10. Repeat the same procedure for each magnification available in your microscope.

6.4 Icons and Functions of the Camera Control Toolbar

Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color / Gray		Video Freeze
	Display Cross Line		Overlay
	Compare stored image with the current live image		Browse Images and Videos in the SD Card
	Settings		Check the LiteView Version

The *Setting* function is relatively more complicated than the other functions. Here are more info about it:

6.4.1 Setting > Measurement


Figure 4 - Measurement Setting Page

<i>Global</i>	<i>Precision</i>	Used for setting digits behind the decimal point for measurement results.
<i>Calibration</i>	<i>Line Width</i>	Used for defining width of the lines for calibration.
	<i>Color</i>	Used for defining color of the lines for calibration.
	<i>EndPoint</i>	<i>Type:</i> Used for defining shape of the endpoints of lines for calibration.
		<i>Null:</i> Means no endpoints.
		<i>Rectangle:</i> Means rectangle type of endpoints. It makes alignment more easily.

Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve
 Left-click the along with the Measurement command mentioned above will unfold the corresponding attribute settings to set the individual property of the Measurement Objects.

6.4.2 Setting > Magnification


Figure 5 - Magnification Settings Page

- Name** The name of the *Magnification*, usually the Magnification of the microscope objective is used as the Magnification name when calibrating, such as *4X*, *10X*, *40X*, *100X*, etc. Besides, other user-defined information could be added into the Magnification name too, for example, microscope model, operator name, etc.
- Resolution** Pixels per meter. Devices like microscopes have high resolution value.
- Clear All** Click the *Clear All* button will clear the calibrated magnifications and resolutions.
- Delete** Click *Delete* to delete the selected item for specific resolution.
- Up** Click *Move Up* to move up the selected Magnification.
- Down** Click *Move Down* to move down the selected Magnification.

6.4.3 Setting > Image Format


Figure 6 - Image Format Settings Page

- Image Format** *JPEG*: The extension of JPEG file can get very high compression rate and display very rich and vivid images by removing redundant images and color data. In other words, it can get better image quality with a minimum least disk space. If measurement objects are available, the measurement objects will be burned into the image and the measurement cannot be edited.
TIFF: Tag Image File Format (TIFF) is a flexible bitmap format that is mainly used to store images including photos and artistic images.
- Measurement Object Save Method** *Burn in Mode*: The measurement objects are merged into the current image. User could not edit the measurement objects any more. This mode is not reversible.
Layered Mode: The measurement objects are saved in different layer with current image data in the target file. User could edit the measurement objects in the target file with some software on the PC. This mode is reversible.

6.4.4 Setting > Video


Figure 7a - Video Settings Page - Playback


Figure 7b - Video Settings Page - Video Encode

Fast Forward / Reverse Interval The time interval of the playback of video files.

Video Encode You can choose *H264* or *H265* encoding. *H265* encoding can significantly reduce encoding bandwidth and save storage space under the same encoding quality.

6.4.5 Setting > Storage


Figure 8 - SD Card Settings Page

Storage Device

SD Card: only SD Card is supported as the storage device.

*File System
Format of the
Storage Device*

List the file system format of the current storage device

FAT32: The file system of SD Card is FAT32. The maximum video file size of single file is 4GB.

exFAT: The file system of SD Card is exFAT. The maximum video file size of single file is 4GB.

NTFS: The file system of SD Card is NTFS. The maximum video file size of single file is 4GB. Use PC to format the SD Card and switch between FAT32, exFat and NTFS.

Unknown Status: SD Card not detected or the file system is not identified.

6.4.6 Setting > Files


Figure 9 - Files Settings Page

Image File Name

Auto: The image files will be saved automatically with the specified prefix.

Manual: User has to specify the file name before image saving.

Video File Name

Auto: The video file will be saved automatically with the specified prefix.

Manual: User has to specify the Video File Name before video recording.

6.4.7 Setting > Language


Figure 10 - Language Selection Settings Page

<i>English</i>	Set language of the whole software into English
<i>Simplified Chinese</i>	Set language of the whole software into Simplified Chinese
<i>Traditional Chinese</i>	Set language of the whole software into Traditional Chinese
<i>Korean</i>	Set language of the whole software into Korean
<i>Thailand</i>	Set language of the whole software into Thai
<i>French</i>	Set language of the whole software into French
<i>German</i>	Set language of the whole software into German
<i>Japanese</i>	Set language of the whole software into Japanese
<i>Italian</i>	Set language of the whole software into Italian
<i>Russian</i>	Set language of the whole software into Russian

6.3.8 Setting > Miscellaneous



Figure 11 - Miscellaneous Settings Page

- | | |
|----------------------------------|---|
| <i>Ruler</i> | Select to display or hide the Ruler in the video window. |
| <i>Measurement</i> | Select to display the Measurement toolbar in the video window, otherwise, the Measurement toolbar will not be displayed. |
| <i>Overlay</i> | Select to support saving graphics Overlay information in fusion mode, or not to save graphics Overlay information in fusion mode. |
| <i>Auto Exposure</i> | The maximum exposure time during auto exposure process could be specified. Setting this item to a lower value could guarantee a faster frame rate during Auto Exposure. |
| <i>ROI Color</i> | Choose the ROI rectangle line color. |
| <i>Camera Parameters Import</i> | Import the <i>Camera Parameters</i> from the SD card to use the previously exported parameters. |
| <i>Camera Parameters Export</i> | Export the <i>Camera Parameters</i> to the SD card to use the previously imported parameters. |
| <i>Reset to factory defaults</i> | <i>Restore</i> camera parameters to its factory status. |


7. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereomicroscopes calibration

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.