ECOVISION Series

INSTRUCTION MANUAL

Model
B-20R
B-20CR
M-100FX
M-100FLed

Ver. 2.0 2020

Summary

1. Warning	3
2. Symbols and conventions	3
3. Safety Information	3
4. Intended use	3
5. Instrument description	4
5.1 B-20R	4
5.2 B-20CR	5
5.3 M-100FX	6
5.4 M-100FLed	7
6. Unpacking	8
7. Assembling	8
7.1 Assembling the microscope	8
7.1.1 Installing the batteries (M-100FLed)	8
8. Use of the microscope	9
8.1 Switching on the microscope	9
8.1.1 B-20R / B-20CR / M-100FX	9
8.1.2 M-100FLed	9
8.2 Light intensity adjustment	9
8.2.1 M-100FX	9
8.2.2 B-20R / B-20CR / M-100FX	9
8.2.3 M-100FLed	9
8.3 Placing the specimen	10
8.3.1 B-20R / M-100FX / M-100FLed	10
8.3.2 B-20CR	10
8.4 Use with batteries	10
8.4.1 M-100FLed	10
8.5 Aperture diaphragm	11
9. Maintenance	12
10. Troubleshooting	13
Equipment disposal	14

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

Standard models


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models


Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. Instrument description


5.1 B-20R


5.2 B-20CR

5.3 M-100FX

5.4 M-100FLed

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

 Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

7.1 Assembling the microscope

- No installation procedure is requested for this microscope series.
- Place the microscope on a sturdy table and connect the power supply.

7.1.1 Installing the batteries (M-100FLed)

1. Using a screwdriver remove the fixing screw of the battery compartment cover. (Fig. 1)
2. Remove the cover (Fig. 2) and insert the batteries (provided), respecting the polarity. (Fig. 3)
3. Put back the battery compartment cover and screw the fixing screw again.

8. Use of the microscope

8.1 Switching on the microscope

8.1.1 B-20R / B-20CR / M-100FX

1. Operate on the main switch ① placed in the rear side of the microscope, moving the selector on "I". (Fig. 4)

8.1.2 M-100FLed

1. Operate on the main switch ② placed in front of the microscope. (Fig. 5)

8.2 Light intensity adjustment

8.2.1 M-100FX

- No light intensity adjustment can be done on this model.

8.2.2 B-20R / B-20CR / M-100FX

1. Turn the light intensity dial ③ placed on the right side of the microscope. (Fig. 6)


8.2.3 M-100FLed

1. Operate on the main switch ② placed in front of the microscope. (Fig. 5)

8.3 Placing the specimen

8.3.1 B-20R / M-100FX / M-100FLed


1. Place the specimen on the stage and secure it using the stage clips ①. (Fig. 7).

8.3.2 B-20CR

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 8)

1. Open the spring arm of the slide holder ② and place the slide from the front on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

8.4 Use with batteries

8.4.1 M-100FLed

- The microscope is provided with an ON-OFF switch ③ (Fig. 9) placed on the bottom side of the microscope. The use of the switch is as follows:
1. **Switch placed in "OFF" position:** use of the microscope is possible only with external power supply. Battery use and recharge is disabled.
 2. **Switch placed in "ON" position:** use of the microscope is possible either with external power supply or with batteries. In this case the battery recharge is enabled.

8.5 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm ring ① (Fig. 10) on the value corresponding to the objective in use. In this case the optimal setting of the condenser is achieved.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 11.


Fig. 10


Fig. 11

9. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

10. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Batteries discharged	Charge the batteries
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor. • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Electric Section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.