

# IM-7



**Inverted Research Microscope**

## IM-7 - Inverted Research Microscope

IM-7 represents the best of what Optika dedicates to the world of research. This model was created to meet all the needs related to research in life science and designed to be complemented by a series of packages dedicated to more advanced individual applications. For all intents and purposes, IM-7 is to be considered as an inverted imaging platform, due to its high expandability and state-of-the-art quality.



# IM-7 - Inverted Research Microscope

Top-level of optical equipment among our product range provides a sharp and clear view in any situation, while top-level mechanical design offers sturdiness and long lifetime.



# IM-7 - Inverted Research Microscope


## IM-7

### **BRIGHTFIELD**

Transmitted brightfield illumination is one of the most commonly used observation method in optical microscopy, and is ideal for fixed, stained specimens or other types of samples having high natural absorption of visible light.

IM-7 is fitted with high-efficiency LED brightfield illuminator, for the best outcome when using this technique.

Capsella middle embryo - Brightfield



### **FLUORESCENCE**

The fluorescence microscopy is the most demanding technique in biology and biomedical sciences, as well as in materials science.

This method is capable to study organic and inorganic samples thanks to primary fluorescence (auto-fluorescence) or secondary (staining and labelling with fluorochromes)

IM-7 is tailored for applications in research, clinical and pharmaceutical diagnostic field.

Fluorescence illuminators available as mercury lamp.

Cotton fibers - UV Fluorescence



## Multiple Observation Methods

### **DIC**

Differential Interference Contrast (DIC) is a microscopy technique that introduces contrast to images of specimens which have little or no contrast when observed using brightfield microscopy. The images produced using DIC have a pseudo 3D-effect, making the technique ideal for many applications.

DIC produces high resolution images with good contrast. It is best for observing unstained samples.

Sphagnum pores g - DIC


### **PHASE CONTRAST**

Phase-contrast microscopy is a particular technique applied in transparent, non-stainable, samples like culture of living cells, microorganisms, lithographic patterns, latex dispersions, fibers, asbestos and subcellular particles.

It reveals many cellular structures that are not visible with a simple brightfield microscope.

Diatoms - Phase contrast

## ② IM-7 - Features



### Revolving nosepiece for DIC and Fluorescence filter turret


The six-position nosepiece has a slot (in each of the six positions) for inserting DIC prisms.

The filter turret can hold up to six fluorescence filterblocks. It is easily extractable, in order to facilitate the operation of inserting or replacing the filterblocks.

### Mechanical stage and universal condenser

The wide 3-layer mechanical stage comes with several interchangeable plates for the use of Petri dishes, flasks and slides. The movement of the stage is controlled by a long tilting handle equipped with a pair of knobs for X/Y axes.

The universal condenser is a 6-position type, designed for brightfield, phase contrast and DIC.



### Illuminator arm

The arm of the transmitted illuminator is backward tilting up to 30 degrees and it allows to use flasks and big bottles.



## Main photo tube and Bertrand lens

The main photo tube located on the binocular head is easily controllable by using its control knob. 3 positions selectable: 100/0, 50/50, 0/100.

For phase contrast centering operations a Bertrand lens is available and it can be easily inserted by means of a dedicated knob.

## Fluorescence equipment



A complete package of accessory dedicated to Fluorescence technique is available as option.



## Fluorescence filterblocks

In addition to the four standard fluorescence B-G-V-UV filterblocks, many others are available upon request to satisfy every kind of need.

## ② IM-7 - Features / Controls





# IM-7 - Features / Controls

- 1 - Main photo tube, on binocular head
- 2 - Handle, for transportation
- 3 - Handle, for transportation
- 4 - Transmitted illuminator switch
- 5 - Fluorescence illuminator switch
- 6 - Magnification changer, 1x-1.5x
- 7 - Right side photo tube
- 8 - Universal condenser, for brightfield, PH and DIC
- 9 - Binocular head, with standard WF10x/22mm eyepieces or WF10x/25mm (optional)
- 10 - Slots for field diaphragm, aperture diaphragm, ND filters

- 11 - Input port for fluorescence illuminator
- 12 - Left side photo tube
- 13 - Side photo tubes control (100/0 ; 20/80 ; 0/100)
- 14 - Main photo tube control (100/0 ; 50/50 ; 0/100)
- 15 - Bertrand lens insertion control
- 16 - Main switch
- 17 - Slot for analyzer
- 18 - Focusing knobs
- 19 - 3-layer mechanical stage
- 20 - Transmitted illuminator brightness control



# IM-7 - Standard Specifications

Part	Description
<b>Optical system:</b>	Infinity corrected.
<b>Head:</b>	Type: Trinocular (Siedentopf); Inclination 45°; Interpupillary distance 47-78 mm; Tube inner diameter 30 mm; Built-in Bertrand lens.
<b>Photo tubes:</b>	Trino Port/Binocular: 100/0 ; 50/50 ; 0/100; Left Side Port/Binocular: 100/0 ; 0/100.; Right Side Port/Binocular: 100/0; 20/80.
<b>Eyepieces:</b>	PLAN WF10x/25 mm, high eyepoint and with built-in rubber cups; focusable .
<b>Nosepiece:</b>	Sextuple ball bearings revolving nosepiece, reversed; DIC slots fo all positions.
<b>Objectives:</b>	Infinity corrected; 60 mm parfocal distance; All with anti-fungus treatment; Selectable according to customer's preferences (see objective table).
<b>Magnification changer:</b>	1x, 1.5x
<b>Specimen stage:</b>	340x230 mm, three-layers mechanical stage; 130x85 mm moving range; Moving mechanism: Rack and pinion (flexible knob); Holder for Petri dish, 160x110 mm; Holder for Terasaki plate (96 well); Holder for 1 slide. Anti-scratch painting.
<b>Focusing:</b>	Coaxial coarse and fine focusing mechanism with limit stop to prevent the contact between objective and specimen; Adjustable tension of coarse focusing knob; 10 mm coarse total travel; Coarse travel per single rotation: 2 mm; Fine travel per single rotation: 0.2 mm; Fine graduations: 100; Fine resolution: 2 µm.
<b>Condenser:</b>	Koehler type, 6-position, focusable, centrable; N.A. 0.55; Iris diaphragm; Phase contrast positions: 10x PH, 20x PH, 40x PH (rings included); DIC positions: DIC1 (10x), DIC2 (20x, 40x, 60x) (prisms not included); 26 mm working distance.
<b>Transmitted illumination (Full Koehler):</b>	X-LED <sup>10</sup> ; 10W LED, high efficiency, 6000K; Brightness control; 65.000-hour lifetime.
<b>Incident illumination:</b>	Upgradable to Incident Fluorescence illumination through the purchase of the components indicated in the fluorescence package table.
<b>Differential Interference Contrast (DIC)</b>	Upgradable to DIC through the purchase of the components indicated in the DIC package table.

# IM-7 - Objectives/Application Packages

## IM-7 is freely configurable in terms of objectives, by choosing among:

Included  Optional

Infinity-corrected Semi-Apochromatic, Long Working Distance objectives, field flatness up to F.N. 25.  
Objectives 20x, 40x, and 60x feature a correction collar in order to compensate for various thicknesses of cover glasses or different containers.

M-1320	IOS LWD U-PLAN F (SEMI-APO) PH 4x/0.13, W.D. 16.5 mm, Cover glass -	<input type="checkbox"/>
M-1321	IOS LWD U-PLAN F (SEMI-APO) PH 10x/0.3, W.D. 7.4 mm, Cover glass 1.2 mm	<input type="checkbox"/>
M-1322	IOS LWD U-PLAN F (SEMI-APO) PH 20x/0.45, W.D. 7.5 - 8.8 mm, Cover glass 0 - 2 mm	<input type="checkbox"/>
M-1323	IOS LWD U-PLAN F (SEMI-APO) PH 40x/0.60, W.D. 3.0 - 4.4 mm, Cover glass 0 - 2 mm	<input type="checkbox"/>
M-1324	IOS LWD U-PLAN F (SEMI-APO) PH 60x/0.70, W.D. 1.8 - 2.6 mm, Cover glass 0.1 - 1.3 mm	<input type="checkbox"/>

## HBO Fluorescence package:

Included  Optional

M-1330	EPI Fluorescence internal attachment	<input type="checkbox"/>
M-151.1	OSRAM 100W HBO high pressure mercury bulb	<input type="checkbox"/>
M-1332	HBO Lamp house	<input type="checkbox"/>
PS-HBO	Optika 100W HBO power supply	<input type="checkbox"/>
M-1334	6-position fluorescence filterbox turret	<input type="checkbox"/>
M-1335	UV protector orange shield	<input type="checkbox"/>
M-1336	B filterblock, filters included	<input type="checkbox"/>
M-1337	G filterblock, filters included	<input type="checkbox"/>
M-1338	V filterblock, filters included	<input type="checkbox"/>
M-1339	UV filterblock, filters included	<input type="checkbox"/>
M-1340	Aperture diaphragm slider	<input type="checkbox"/>
M-1341	Field diaphragm slider	<input type="checkbox"/>
M-1342	Slider with neutral density filter for HBO illumination	<input type="checkbox"/>
M-1343	Empty fluorescence filterblock	<input type="checkbox"/>

Name	Excitation filter (nm)	Dichroic cut-off mirror (nm)	Emission filter (nm)
<b>B (Blue)</b>	460 - 490	500	510LP
<b>G (Green)</b>	510 - 550	570	590LP
<b>V (Violet)</b>	400 - 410	455	455LP
<b>UV (Ultraviolet)</b>	330 - 385	400	420LP

## DIC - Differential Interference Contrast package:

Included  Optional

M-1350	DIC Prism for 10x	<input type="checkbox"/>
M-1351	DIC Prism for 20x	<input type="checkbox"/>
M-1352	DIC Prism for 40x	<input type="checkbox"/>
M-1353	DIC prism for 60x	<input type="checkbox"/>
M-1354	DIC 1 prism 10x for condenser	<input type="checkbox"/>
M-1355	DIC 2 prism 20x-40x-60x for condenser	<input type="checkbox"/>
M-1356	Slider with rotating analyzer	<input type="checkbox"/>

# IM-7 - Accessories

## Eyecups & Eyepieces

- M-1360 PL10x/25, high eyepoint, focusable, rubber cup  
 M-1361 PL10x/25, micrometric eyepiece, high eyepoint, focusable, rubber cup

## Objectives

- M-1320 IOS U-PLAN F (SEMI-APO) PH 4x/0.13, W.D. 16.5 mm, Cover glass -  
 M-1321 IOS U-PLAN F (SEMI-APO) PH 10x/0.3, W.D. 7.4 mm, Cover glass 1.2 mm  
 M-1322 IOS U-PLAN F (SEMI-APO) PH 20x/0.45, W.D. 7.5 - 8.8 mm, Cover glass 0 - 2 mm  
 M-1323 IOS U-PLAN F (SEMI-APO) PH 40x/0.60, W.D. 3.0 - 4.4 mm, Cover glass 0 - 2 mm  
 M-1324 IOS U-PLAN F (SEMI-APO) PH 60x/0.70, W.D. 1.8 - 2.6 mm, Cover glass 0.1 - 1.3 mm

## Fluorescence package

- M-1330 EPI Fluorescence internal attachment  
 M-151.1 OSRAM 100W HBO high pressure mercury bulb  
 M-1332 HBO lamp house  
 PS-HBO 100W HBO power supply)  
 M-1334 6-position fluorescence filterbox turret  
 M-1335 UV protector orange shield  
 M-1336 B filterblock, filters included  
 M-1337 G filterblock, filters included  
 M-1338 V filterblock, filters included  
 M-1339 UV filterblock, filters included  
 M-1340 Aperture diaphragm slider  
 M-1341 Field diaphragm slider

## DIC package

- M-1350 DIC Prism for 10x  
 M-1351 DIC Prism for 20x  
 M-1352 DIC Prism for 40x  
 M-1353 DIC Prism for 60x  
 M-1354 DIC 1 prism 10x for condenser  
 M-1355 DIC 2 prism 20x-40x-60x for condenser  
 M-1356 Slider with rotating analyzer

## Camera Adapters

- M-620 0.35x focusable C-Mount adapter (for main photo tube only)  
 M-620.1 0.5x focusable C-Mount adapter (for main photo tube only)  
 M-620.2 0.65x focusable C-Mount adapter (for main photo tube only)  
 M-620.3 1x focusable C-Mount adapter (for main photo tube only)  
 M-1365 0.5x focusable C-Mount adapter (for left/right side photo tube only)  
 M-1366 1x focusable C-Mount adapter (for left/right side photo tube only)

## Miscellaneous

- M-1370 Color temperature filter, 38mm  
 M-1371 Green filter, 38mm  
 M-1372 Yellow filter, 38mm  
 M-1373 Frosted filter, 38mm  
 15104 Cleaning kit  
 DC-005 TNT dust cover, extra large, 820(l)x550(h) mm  
 M-005 Micrometric slide, 26x76mm, with 2 scales (1mm/100 & 10mm/100)  
 VP-IM7 IQ/OQ/PQ manual for IM-7  
 AB-040 Antibacterial surface treatment, only for newly purchased microscope

ÖPTIKA



Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
 Tel: (877) 877-7274 | Fax: (516) 801-2046  
 Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)



### How to connect the cameras to our microscopes.

Please refer to the Adapter reference list on Digital section.

### 15104 - Cleaning kit

It cleans glass quickly and effectively, without leaving residue or odor. Ideal for precision lens or prism cleaning.



v 6.7 - OPTIKA reserves the right to make corrections, modifications, enhancements, improvements and other changes to its products at any time without notice.

## Headquarters and Manufacturing Facilities

**OPTIKA® S.r.l.** Via Rigla, 30 - 24010 Ponteranica (BG) - ITALY - Tel.: +39 035.571.392 - [info@optikamicroscopes.com](mailto:info@optikamicroscopes.com)

## Optika Sales branches

**OPTIKA® Spain** [spain@optikamicroscopes.com](mailto:spain@optikamicroscopes.com)  
**OPTIKA® China** [china@optikamicroscopes.com](mailto:china@optikamicroscopes.com)  
**OPTIKA® India** [india@optikamicroscopes.com](mailto:india@optikamicroscopes.com)

**OPTIKA® USA** [usa@optikamicroscopes.com](mailto:usa@optikamicroscopes.com)  
**OPTIKA® Central America** [camerica@optikamicroscopes.com](mailto:camerica@optikamicroscopes.com)