# OPTIKA


M I C R O S C O P E S  
I T A L Y

## XDS - 1R

## DESCRIZIONE

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# INVERTED MICROSCOPE

## OPERATION INSTRUCTION

The model XDS-1R Inverted Microscope has the feature that can be observed in culture bottle or culture utensil, it is furnished with the long working distance plan achromatic objective, long or ultra-long working distance condenser, wide field eyepiece and phase contrast attachment. It can be observed transparent and unfinished living specimen. It is advanced instrument suitable for research of living specimen structure, liquid, deposits and so on. It also can be used in research institutes, universities, medical treatment, agriculture and animal husbandry etc.

## I . SPECIFICATION

### 1.Eyepiece

Type	Magnification	Foucs (mm)	Field (mm)	Remark
Wide field eyepiece	10×	25	Φ 18	
Plan eyepiece	16×	15.6	Φ 11	select

### 2.Objective

Type	Magnification	N.A	W.D (mm)	Cover glass thickness (mm)	Remark
Long working distance plan achromatic objective	10×	0.25	8.1	—	
	25×	0.4	4.8	1.2	
	40×	0.60	3.3	1.2	
Long working distance plan phase contrast objective	10×	0.25	8.1	—	
	25×	0.4	4.8	1.2	Select
	40×	0.60	3.3	1.2	select

### 3. Total Magnification

Total Magnification	Objective Magnification	10×	25×	40×
Eyepiece	10×	100×	250×	400×
	16×	160×	400×	640×

#### 4. Condenser

- (1) long working distance condenser (with phase contrast): NA=0.4,  
working distance : 28mm;
- (2) ultra-long working distance(select): working distance: 70mm.

#### 5. Stage

size: 200mm × 152mm;  
movable range: 30mm × 70 mm.

#### 6. Coaxial coarse/fine focusing system with limit knob and tensional adjustment ring


minimum division of fine focusing: 0.002 mm.

#### 7. Interpupillary distance : 53~75mm.

#### 8. Illumination system: 6V,20W halogen bulb, light intensity adjustable.


#### 9. Power supply: 220V(50HZ) or 110V(60HZ).

## II. COMPONENTS

**XDS-1R (Fig.1)**

- 1.lamp box    2.lamp box holding screw    3.filter base    4.annular  
diaphragm adjustment screw    5.annular diaphragm slab    6.eyepieces  
7.trinocular    8.tensional adjustment ring    9.coarse focusing knob  
10.fine focusing knob    11.objective    12.stage    13.condenser adjustment  
screw    14.aperture diaphragm swivel    15.lamp holder adjustment handle

### XDS-1R (Fig.2)

16.collector adjustment handle 17.condenser holding screw 18.condenser up/down knob 19. lengthwise movement knob 20. cross movement knob 21.power switch 22. light control 23. limit knob 24.tube holding screw 25. observation or photograph lever 26.field diaphragm adjustment screw

### III. INSTRUMENT OPERATION

1. Insert the condenser into the lift bracket, plug into the outlet, hold condenser by condenser holding screw 17, then push power switch 21 into "I", power supply is switch on. (Fig.2)
2. Put the specimen on the center of stage, turn  $10\times$  objective to working position, then adjust coarse / fine focusing system.
3. Adjust interpupillary distance and diopter.
4. Adjust the position of condenser, light control 22 and aperture diaphragm swivel 14 to gain satisfactory illumination state.(Fig.1,2)
5. Little adjust fine focusing knob 10 when use different magnification objective.(Fig.1)

### IV. COMPONENTS DETAILED OPERATION :

#### 1. Adjust of interpupillary distance (Fig.3)

Put specimen on stage and ring specimen into exact focus. Adjust interpupillary distance of binocular until right-left field of view can be composed one.


Fig.3

#### 2. Adjustment of diopter(Fig.3)

Put specimen on stage, turn  $40\times$  objective to working position. Firstly, observe at right tube with right eye, adjust coarse / fine focusing knob to image clearly. Secondly observe at left tube with left eye, adjust diopter ring 1 to image clearly.


Fig.4

#### 3. Coarse /. fine focusing

The instrument used coaxial coarse/fine focusing mechanism. Tensional adjustment ring 3 used for adjusting tension of coarse focusing knob 2 to prevent stage from naturally sliding down. Limit knob 1 prevents accidental

specimen/objective contact. (Fig.4)

#### 4. Stage

Two different bedplate (1 or 4) can be put on stage, used for different culture and culture utensil, lengthwise / cross is very convenient for adjustment.(fig.5)


Fig.5

#### 5. Adjustment of lamp position

Turn condenser up/down knob 18 to move long working distance condenser to reticle position, then put a piece of white paper on base of filter, loosen lamp holder holding screw 2 (Fig.8), move lamp holder adjustment handle 15 and collection adjustment handle 16 to the imaging of filament be on the white paper, then tighten the lamp holder holding screw, fixup the position of lamp. (fig.1, 2)

#### 6. Adjustment of field diaphragm

Turn 10×objective to working position, observe with 10× objective, turn coarse / fine focusing knob 9 and 10 to make the image of specimen clearly. Turn field diaphragm adjustment handle 26 to diminish field diaphragm, turn condenser up/down knob 18 to gain the image of field diaphragm, then adjust condenser adjustment screw 13 to center field diaphragm and eyepiece diaphragm, turn field diaphragm adjustment handle to make field diaphragm less than eyepiece diaphragm.(Fig.1,2)

#### 7. Adjustment of condenser aperture diaphragm

Turn aperture diaphragm swivel 14 to suit the numerical aperture of condenser with the numerical aperture of objective, in the way, can gain

opportune contrast image and satisfactory illumination.(Fig.1)

#### 8. Power switch and light control

Push power switch 21 into “ I ”, power is switched on, then adjust light control 22 to observe the image of specimen comfortably.(Fig.2)

**NOTICE:** Don't put the light control on the light highest position long, avoid depressing life of lamp.

#### 9. Ultra long working distance condenser (select)

Loosen condenser holding screw before using ultra long working distance condenser, take out long working distance condenser, then install ultra long working distance condenser and hold it. Turn condenser up/down knob to make the facula of specimen focusing into light point.

#### 10. Phase contrast attachment

- (1) Adjust field diaphragm according to step 6, make long working distance condenser move to corresponding position.
- (2) Insert annular diaphragm slab 5 in condenser, turn aperture diaphragm swivel 14 to make aperture diaphragm maximal.(Fig.1)
- (3) Turn corresponding multiple phase contrast objective to working position.
- (4) Take out one eyepiece, insert center telescope in eyepiece tube and adjust it, make the image of annular phase plane and annular diaphragm clear.(Fig.6)
- (5) Adjust annular diaphragm adjustment screw 4, make annular diaphragm light halo and annular phase plane dark halo superposition.(Fig.7)


Fig.6

(6) Take out center telescope and insert eyepiece to observe phase contrast.

Note: It is a necessary process to revise annular diaphragm observing with phase contrast.

11. Exchange the lamp and fuse

Exchanging the lamp and fuse according to the following steps (Fig.8)

(1) Switch off power supply and pull out plug.

(2) loosen lamp holder holding screw 2 and lamp adjustment handle 1, remove lamp holder 3 from lamp box.

(3) Pull out the old lamp from lamp holder 1, then insert the new lamp into lamp holder, and clean the new lamp with alcohol.

(4) Reinstall lamp holder into lamp box, switch on power supply again, move lamp adjustment handle, observe the uniformity of the illumination, then hold the lamp holder by lamp holder


Fig.7


Fig.8


Fig.9

holding screw.

(5) Loosen screw of fuse 1, put out the bad fuse, mount the new fuse, and tighten screw of fuse.

the specification of fuse:  $\Phi 5$ , 0.5A(110V) or 0.25A (220V). (Fig.9)

## IV.MAINTENANCE

### 1. Sweep the lens

Sweep the lens by lens tissue or soft fabric immersed with mixed liquid of alcohol/ether or diethyl benzene.

### 2. Clean the painted parts

The dust on the painted parts can be removed by gauze, for the grease spots, the gauze immersed slightly with aviation gasoline is recommended. Do not use organic solvents such as alcohol, ether or other thinner etc, for cleaning the pointed parts or plastic components.

### 3. Avoid disassembling the microscope

Being a precise instrument, do not disassemble the microscope casually that may cause serious damage to its performance.

### 4. Being not used

Cover the microscope with polymethyl methacrylate or polyethylene and places where there is dry and mouldless. Suggest that storage all objectives and eyepieces in closed container with drying agent.

**Note: Buying lamp come to our factory**

**OPTIKA S.R.L.**  
Via Rigla 30, Ponteranica (BG) - ITALY  
Tel.: ++39 035 571392 (6 linee) Telefax: ++ 39 035 571435

**MAD Iberica Aparatos Cientificos**  
c/. Puig i Pidemunt, nº 28 1º 2<sup>a</sup> - (Pol. Ind. Plà d'en Boet) 08302 MATARÓ  
(Barcelona) España Tel: +34 937.586.245 +34 937.414.529

**New York Microscope Company Inc**  
100 Lauman Lane, Suite A, Hicksville, New York 11801, USA  
Tel.: 877.877.7274 - Fax: 516.801.2046  
[www.microscopeinternational.com](http://www.microscopeinternational.com) - [info@nscopes.com](mailto:info@nscopes.com)