# XDS-2 + M-795

1.0 COMPONENTS	page 6
2.0 INSTALLATION	page 8
3.0 ADJUSTMENT	page 13
4.0 ADJUSTING THE MICROSCOPE	page 14
5.0 PHASE CONTRAST OPERATION	page 18
6.0 LAMP MAINTENACE	page 19
7.0 MICROPHOTOGRAPHY AND VIDEO	page 21
8.0 TECHNICAL SPECIFICATION	page 23
9.0 TROUBLESHOOTING	page 24
10.0 RECOVERY AND RECYCLING	page 26

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this ma-nual.

#### Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

#### **Electrical safety**

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

#### Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

The instructions on this symbol to avoid possible severe personal injuries.

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

Possibility of electric shock.


Attention: high temperature surfaces. Avoid direct contact.


Technical notes or usage tips.

#### **SAFETY NOTE**


- 1 Do not store the instrument under direct sunlight, high temperature or high humidity, nor in dusty environments. Make sure the stage is plane, horizontal and stable.
- When moving the microscope, please hold the instrument with one hand on the lower side of the eyepiece tube (1), and the other hand on the illumination bracket (2) (Figure 1).
- 3 If water or bacterial cultures are spilled on the stage, objective or viewing tube, pull out the power cord at once, and wipe the microscope clean. Failing to do so may result in damage to the instrument.

**7 XDS-2FL:** Be sure the lamp is cool before you cover the microscope with the dust-cover after use.

**9 XDS-2FL:** The lamp light contains ultraviolet radiation that is harmful for the eyes and skin. Always look at the lamp light through the orange shield.

#### **MAINTENANCE**

- 1 Use a lint-free cloth to wipe the glass parts. In order to remove fingerprints and oil stains, slightly dampen the cloth with xylene or with a 3:7 mixture of ethanol and ether.
- Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- **2** Do not use organic solvents to wipe the non-optical elements. Please use a neutral detergent only.
- **3** During use, if liquids fall on the microscopes, cut off the power at once, and wipe up the moisture.
- **4** Do not disassemble the microscope. Disassembling might degrade the performance of the microscope, and will invalidate the warranty.
- When the objectives are not mounted, please cover the sockets with the dust caps, in order to prevent dust and liquids from entering the microscope.

Figure 1

Figure 2

When not in use, remember to cover the microscope with the dust hood. Please let the lamp cool down before covering the microscope.

#### **SAFETY SYMBOLS**

Symbol	Meaning
	The surface may be hot, do not touch.
$\overline{\bigvee}$	Please read the instruction carefully, before use. Improper operation will result in injuries or malfunction.
_	Main switch – ON position.
0	Main switch – OFF position.

#### Connect the mains plug into the socket at the base

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.


The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.


#### AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.


#### 2.1 INSTALLATION DIAGRAM

The following figure shows the installation sequence of the components.

The number in the figure shows the installation steps.

- ▶ Before installing, be sure every component is clean. Take care not to scratch parts or surfaces.
- ► Store the supplied hexagon wrench in a safe place. When changing the components, you will need it again.

#### 2.2 INSTALLATION STEPS

#### 2-2-1 Installing the condenser illumination unit (fig. 3)

- 1 Insert the condenser illumination unit (3) into the bracket, according to figure 3.
- 2 Turn the condenser illumination unit clockwise about 90°, with the "AS" mark of filter holder (3) facing forwards. Align the screw of the condenser illumination unit and the hole of the holder, then screw the bolt in the hole with the supplied hexagon spanner.
- **3** Insert the connector plug (4) into the connector jack (5).

#### 2-2-2 Installing the LED housing (figure 4)

Align the connector plugs (1) and the LED housing pins (2) Push the LED housing gently into the illumination unit.

#### 2-2-3 Installing the objectives (fig. 5 & fig. 6)

- 1 Turning the coarse focusing knob (1) till the nosepiece reaches to its lowest position.
- For a safe transport, the nosepiece is placed in the lowest position and the tension adjustment collar (2) is adjusted to the appropriate tension when the microscope leaves the factory.
- 2 Screw the lowest magnification objective on to the turret from the right side, then turn the turret clockwise. Mount the other objectives in the same way, following the sequence from low to high.
- ▶ **Note**: the objectives can also be installed through the stage opening.
- ► Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.
- ➤ To prevent dust and contamination from entering the microscope, cover all the unused holes with dust caps (3).
- ► When operating, use the low magnification objective (4x or 10X) to search and focus the specimen, then switch to higher magnifications.
- ▶ When switching between objectives, slowly turn the nosepiece until it clicks. The click means that the objective is in the right position, in the center of the light path.

Figure 3

Figure 4

Figure 5

Figure 6

#### 2-2-4 Installing the stage extension and the mechanical stage (fig. 7)

- ► The stage extension can be installed on either side of the stage to enlarge the working surface. The mechanical stage must be installed on the side opposite the extension.
- ► For right-handed operators, the mechanical stage is normally installed on the right side.
- 1. Installing the stage extension: First, screw the bolts (1) on to the extension, then mount the extension from below the stage. Screw it to a firm fit.
- 2. Installing the mechanical stage: As for the extension, the mechanical stage is fixed by two bolts under the stage.

Figure 7

#### 2-2-5 Installing the stage insert (fig. 8)

- **1.** When using the glass stage (1), make sure that the stage is horizontal.
- **2.** Install the stage insert in the stage opening.
- ► Turn the disk until the V-groove faces the user. This simplifies objective identification.


Figure 8

#### 2-2-6 Installing the eyepieces (fig. 9)

- **1.** Remove the cap of the eyepiece tubes (1).
- 2. Insert the eyepieces into the tubes.


Figure 9

#### 2-2-7 Installing the color filters (fig. 10)

- ▶ Be sure the color filter has cooled down completely before you change them. Remove the filter holder (1), then install the color filters (2) you need.
- Mount the color filter flat as shown in (3), verifying that they are not tilted.
- ► If the color filter is tilted or otherwise out of place (4), it may fall.
- The color filters can be stacked in the holder. This allows to install as many filters as needed, as long as the whole thickness is less than 11 mm.

#### 2-2-8 Connecting the power cord (fig.11, 12 & 13)

- ▶ Do not tension or otherwise stress the power cord. If bent, the wires may easily break.
- **1.** Turn the main switch (1) to "O"(off) before connecting the power cord.
- **2.** Insert the plug (2) into the power jack (3) of the microscope.
- **3.** Plug the power cord (4) into the mains socket. Check for a safe connection.
- 4. Insert the low voltage connector plug (9) into the socket (8).
- ► Please use the supplied power cord. If lost or damaged, please refer to qualified service.

Figure 11

#### 2-2-9 Replacing the fuse (fig.11 & 12)

Before replacing the fuse, turn the main switch (1) to "O" (off) and unplug the power cord. Rotate the fuse (6) support out of the holder (7) using a straight screwdriver. Insert a new fuse in the support, then rotate the support back into the holder.

Fuse rating: 250V, 500mA.

Figure 12

Figure 13


# 2-2-10 XDS-2FL: Insertion of the mercury lamp housing (fig. 13 & 14)

Put the lamp housing behind the microscope and fully insert it in the lodging (1).
With the provided Allen wrench, tighten the two screws on the top of the lodging, (2) and (3).

#### 2-2-11 XDS-2FL: Insertion of filter holder (fig. 15)

You can insert the filter holder in front of the mercury lamp,

in the appropriate slot.

Put the protection shield on the illumination tube just above the stage and lock it in position.

# 2-2-13 XDS-2FL: Connection of the fluorescence power supply (fig. 17)

On the back of the NFP-1 power supply, plug in the 3-poles connector (that comes from the lamp housing) and lock it by screwing the ring nut (1). Then plug the line voltage in (2).

Pay attention to the right position of the red line-voltage selectors on the back of NFP-1: they must match your line voltage specification: 50/60Hz; 115/230Vac

Figure 14

Figure 15


Figure 16

Figure 17

## 4.0 ADJUSTING THE MICROSCOPE

#### **4-1 MICROSCOPE BASE**

#### 4-1-1 Turning on the LED (fig. 18)

Connect the power, turn on the main switch (1) (fig. 18).

#### 4-1-2 Adjusting the brightness (fig. 19)

Turn the brightness adjustment knob clockwise to increase the brightness, and counteclockwise to decrease it.

Figure 18

#### 4-1-3 Adjusting the tension (fig. 20)

► The coarse focusing knob (2) is preadjusted to a tight tension upon leaving the factory.

If the nosepiece drops down by itself, or the specimen defocuses while adjusting the fine focus knob (3), the coarse focus knob is too loose.

Turning the tension adjustment collar (1) in the indicated direction tightens the coarse focus tension (2). Rotate in the opposite direction to decrease the tension.

Figure 19


Figure 20

#### 4-2 STAGE

#### 4-2-1 Setting the specimen (fig. 21 & 22)

- For the best image quality, use flasks, Petri dishes and slides with a 1.2 mm thickness.
- Using Φ35mm culture dishes: You can place a Φ 35 mm culture dish on the stage by using the standard center board (1) of the stage.

- **1.** When observing 96-well or 24-well microtitration plates, please fasten them directly with the stage clip (2).
- **2.** When fastening other plates, please use the adaptor plates supplied with the mechanical stage:
- Terasaki bracket (3) for Terasaki boards
- Culture dish bracket (4) for Φ 35 mm culture dishes
- Object slide bracket (5) for object slides and Φ 54 mm culture dishes

Figure 21

Figure 22

**3.** Turning the X and Y knobs (6,7), move the specimen to the required position. (Movement Range: 120 (width) × 78 (length) mm).

#### 4-2-2 Moving the specimen

Move the specimen to the desired position by freehand or by turning the knobs of the mechanical stage.

► When switching objectives, take care not to touch the adaptor plates with the objectives, as their weight may damage the front lens.

#### **4-3 VIEWING TUBE**

#### 4-3-1 Dioptric adjustment (fig. 23)

- 1. Look into the right eyepiece with your right eye only, and focus on the specimen.
- 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring (1) to compensate.
- ► The adjustment range is ±5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction (fig.23).

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

► The graduation on the interpupillary distance indicator (3), pointed by the spot "." (2) on the eyepiece holder, shows the distance between the operator's eyes.

The range of the interpupillary distance is 48 ~ 75mm.

### 4-3-3 Selecting the light path (fig. 25)

 Pull the light path selector lever (1) sideways using your thumb, selecting the light path you need.

LIGHT PATH SELECTOR LEVER	BRIGHTNESS	APPLICATION
IN	100% used for binocular observation	Binocular observation
OUT	20% used for binocular observation, and 80% used for video or photography	Binocular observation television, and micrography or video can be operated simultaneously

Figure 23

Figure 24

Figure 25

#### 4-4 ILLUMINATION UNIT

#### 4-4-1 Using color filters (fig. 26)

- Selecting the appropriate color filtersaccording your need.
- You can stack a group of color filters in the filter holder, if you ensure that they are level and that the whole thickness is less than 11mm.

COLOR FILTER	USE
IF550	Single contrast color filter (green) (used for phase contrast microscopy)
LBD	Color temperature compensation color filter (Blue) (used for bright field observation and microphotography)


Figure 26


Figure 27

#### 4-4-2 Using the aperture diaphragm (fig. 27)

- ▶ When in bright field observation, the aperture diaphragm controls the numerical aperture of the illumination system. When the numerical aperture of the objective and the aperture of the illumination system match, the highest resolution is achieved.
  - To identify the aperture diaphragm, the eyepiece can be removed or the centering telescope can be used. The aperture diaphragm looks as in fig. 27. The aperture can be changed by moving the aperture adjustment lever. ((1) is the image of the aperture diaphragm, (2) is the edge of the objective).
  - Generally, when observing a fully chromatic specimen, you need to set the size of the condenser to 70-80% of the aperture of the objective. When observing unstained samples (e.g. bacteria), start from 70% and slowly turn the aperture diaphragm lever clockwise.

#### 5-1 IDENTIFYING THE COMPONENTS

#### 5-1-1 Phase contrast objectives (fig. 28)

- The magnifications of the optional phase contrast objectives are: 10X, 20X. Such objectives are marked "PH".
- Mount the objectives on the turret using the same procedures as the standard objectives (2-2-3).

- Adjustable phase slider.
- The light ring is pre-centered when the microscope leaves the factory. It should therefore need no further adjustment. If a recentering is needed, it can be performed via the two side bolts.
- The 10X/20X light ring (1) must be used with the 10X, 20X phase contrast objectives, while the opening (2) is used for bright field.

Figure 28

Figure 29

Figure 30

#### **5-2 SETUP AND USE**

#### 5-2-1 Installing the phase contrast slider (fig. 30)

- **1.** Insert the slider (1) into the illumination system, printed face up.
- 2. Pull the slider into the desired position, to the click stop.
- 3. When in phase contrast observation, keep the aperture diaphragm adjustment lever (2) on the "O" (open) position.

### 6-1 Aligning the lamp (fig. 31)

Move the filter block to the desired position.

Move the lamp using the two lamp centring knobs on the side (1) until no dark zones are visible on the specimen you are viewing.

Then adjust the lamp focus knob (2) until you can obtain an even illumination of the whole field of view. You can optimize the illumination by adjusting the tilt angle of the rear mirror of the lamp (3).



Figure 31

#### 6-2 Selecting the filter set (fig. 32)

Switch between Green and Blue excitation filter sets by pulling the bar (1) under the objectifs. "BF" position allows to view in normal brightfield mode.

#### 6-3 Maintenance and warnings

- Always turn on the lamp with the current adjusted to maximum (5A) and then wait for a minimum of 10 minutes before using the microscope.
- If the maximum current goes down below 4.5 A the lamp needs to be changed.
- After the lamp has been switched off it takes about 30 minutes before it is cool enough to be ready to be used again.
- Be sure the lamp is cool before you cover the microscope with the dust-cover after use.
- The lamp has a lifetime of 500 hours. In order to avoid explosive failure, replace the lamp after 500 hours, even if it seems to work correctly.
- The lamp light contains ultraviolet radiation that is harmful for the eyes and skin. Always look at the lamp light through the orange shield or through the dark alignment window on the filter block.

#### 6-3-1 Mercury lamp insertion and replacement (fig.33)

Disconnect any cable from the epiilluminator. Loosen the lock-screw (1) on the lamp housing. The lamp holder will come out with the door. Remove the old lamp or the plasticplaceholder. Use protection glasses as, if handled carelessly, there is a risk that the lamp may shatter and eject glass splinters.

Thus, without touching the glass part of the lamp, take out the fluorescence lamp from its box, taking care not to exert stress on any part of the lamp. Put one end inside the flexible metal holder (observe that it fits only in one direction) and then in the fixed one.

Tighten the screws just enough to make it stay in place. If you tighten them too much the lamp may shatter. Put the lamp door and lamp assembly back in place and connect the cable on the side of the epiilluminator to the power supply.


Figure 33

Connect the mains cable of the power supply. Turn on the power supply using the power switch. You should now wait for at least 10 minutes before aligning the lamp.

On the life time display can be read for how many hours the lamp has been in use. Having replaced the lamp, press the reset button to reset the life time count.

## 7.0 MICROPHOTOGRAPHY AND VIDEO

#### 7-1 VIDEO PORT (VIDEO)

#### 7-1-1 Selecting the light path (fig. 34)

**1.** To activate the video port, pull out the light path selector lever.

For observation of dark specimens, first focus, then pull out the lever.

Figure 34

#### 7-1-2 Installing the video adapter (fig.35)

- 1. Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
- **2.** Set up the video adapter (3) on the camera according to its instructions.
- **3.** Install the adapter into the tri-through port, and screw down the bolt (1).

#### 7-1-3 Focus (fig. 35)

During a binocular observation at 20% brightness, look at the image on the video imaging system, refocusing the adapter if necessary (4).


Figure 35

### 7.0 MICROPHOTOGRAPHY AND VIDEO

#### 7-2 MICROPHOTOGRAPHY

#### 7-2-1 Selecting the light path

Set up the light path as in 7-1-1 and in 4-3-3.

#### 7-2-2 Installing the photography adapter (fig. 36)

- **1.** Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
- 2. Install the photography adapter (3) into the tri-through port according to its instructions, and screw down the locking bolt (1).
- 3. Inserted the camera head (4) (if any) into the adapter.
- **4.** Plug the adapter+head assembly into the photo tube, then screw down the locking bolts (1).
- Connect the camera+head assembly to the camera according to the camera instructions.
- For the photography of dark specimens, obscure the eyepieces and the viewfinder with a dark cloth in order to reduce stray light.
- The camera magnification can be calculated as objective magnification × camera + lens magnification.
- ► When shooting with a SLR, the mirror movement may cause camera movement. Please lift the mirror, use long exposure times and use an extension cord.


Figure 36

#### 7-2-3 Focus

During a binocular observation at 20% brightness, look at the image on the video imaging system, refocusing the adapter if necessary.

#### 7-2-4 Adjusting the color temperature

When using a film camera, please use sunlight film.

- 1. Mount the LBD color temperature correction (blue) filter to the color filter bracket.
- 2. Turn the brightness adjustment knob to the maximum position.

### **8-1 MAIN SPECIFICATIONS**

Illumination	Light source type X-LED <sup>8</sup> with white 8W LED; light intensity control using a knob on front side of the frame.  Color temperature: 6300K  LED average life time approx. 50.000h  Voltage: 110/240Vac, 50/60Hz, 1A; Fuse: T500mA 250V  Max power required: 13W	
<b>Observation Modes</b>	Brightfield, phase contrast.	
Focusing	Coaxial coarse and fine focusing mechanism (graduated, 0.002mm) with upper stop, to prevent the contact between objective and specimen.  Adjustable tension of coarse focusing knob.	
Stage	Fixed stage, dimensions 250x230 mm.  Mechanical stage mountable on the right side of the stage, X-Y translation range 119x70 mm, with metallic interchangeable inserts for slides, Petri dishes, Terasaki, multi-Well plates, etc.  Pair of side extensions to expand the surface of the stage.  Glass stage insert with hole for small dimension specimens.	
Nosepiece	Quintuple revolving nosepiece, rotation on ball bearings.	
Head	Trinocular observation head, inclined 30° and rotatable 360°. Diopter adjustment on left eyepiece. Interpupillary adjustment 55-75 mm.	
Eyepieces	Wide field eyepieces EWF10X/22 with field number 22.	
Objectives	Infinity corrected optical system IOS (Infinity Optical System). Plan-achromatic LWD objectives infinity corrected, for thickness 1,2 mm, made by following objectives: -) Plan-achromatic IOS LWD 4X, N.A. 0.10, W.D. 18,0 mm -) Plan-achromatic IOS LWD 10XPh, N.A. 0.25, W.D. 10,0 mm -) Plan-achromatic IOS LWD 20XPh, N.A. 0.40, W.D. 5,1 mm -) Plan-achromatic IOS LWD 40X, N.A. 0,65, W.D. 2,6 mm All objectives are treated with an anti-fungus treatment.	
Condenser	LWD condenser, N.A. 0.30, working distance 72 mm. The condenser can be removed to extend the working distance up to 150 mm. Precentered slider with phase ring 10X/20X.	
Dimensions	HEIGHT: 473 mm WIDTH: 250 mm DEPTH: 515 mm WEIGHT: 9 kg	
Accessories	Instruction manual and dust cover included	

### **8-2 OBJECTIVE SPECIFICATIONS**

TYPE	MAGNIFI- CATION	NUMERICAL APERTURE (N.A.)	WORKING DISTANCE (mm)	CONJUGA- TE DISTANCE (mm)	FOCUS DISTANCE (mm)	COVER SLIP THICKNESS
Infinite Long Working Distance Plan Achro-	4X	0.1	25.2	00	45	-
matic Objective	40X	0.6	3.2	00	45	1.2mm
Infinite Long Working Distance Plan Phase Contrast Objective	10X	0.25	11	00	45	0.17
	20X	0.4	6	00	45	0.17

Under certain conditions, some no-fault factors will bring a reversible influence to the instrument's performance. If this happens, please take proper measures according to the follow table. If you can't solve the trouble as indicated, please contact the sales department of our Company.

PROBLEM	REASON	SOLUTION
I. Optical section		
1. The illumination is open, but the field of view is dark.	The plug of the LED holder is not connected to the illumination set	Connect them
	The brightness is too low	Adjust to a proper setting
	Too many colour filters have been stacked	Minimize the number of the filters
2. The edge of the field of view is vignetted or the brightness is asym-	The nosepiece is not in the correct position	Turn the nosepiece to a click stop
metric.	The color filter is partially inserted	Insert the filter to full depth
	The phase contrast slider is not in the proper position	Move the slider to a click stop
3. Dust and stains can be seen in the field of	There are stains and dust on the specimen	Clean the specimen
view.	There are stains and dust on the eyepiece	Clean the eyepiece
4. There is an apparent double image.	The size of the aperture diaphragm is too small	Open the aperture diaphragm
5. Poor image quality: The image is not sharp	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
The contrast is not high The details are not clear The phase contrast is low.	The aperture diaphragm in the view of field is opened too much or too little	Adjust the aperture diaphragm
	The lenses (condenser, objective, eyepieces are culture dish) is dirty	Thoroughly clean all the optical system
	In phase contrast observation, the bottom thickness of the sample is more than 1.2mm	Use a sample holder whose bottom thickness is less than 1.2mm
	A bright field objective is used for phase contrast observation	Switch to a phase contrast objective
	The condenser ring is not aligned with the objective phase ring	Adjust the condenser ring to match the objective phase ring
	The light ring and/or the phase contrast ring is not centered	Adjust the bolts to center them
	The objective used is not compatible with the phase ring	Please use a compatible objective
	The phase contrast depends on the sample position	The sample holder is not flat. Move the sample around until a compatible area is found.

# 9.0 TROUBLESHOOTING

6. One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical section		
1.The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
2.The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section		
1. The LED doesn't turn on.	No power supply	Check the power cord connection
2. The brightness is not enought	The brightness adjustment is low	Adjust the brightness
3. The light blinks	The power cord is poorly connected	Check the power cord
IV. Viewing tube assemb	ly	
The field of view of the two eyes is different.	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptric correction is not right	Adjust the dioptric correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and	d video	
1. The image is unfocused.	Incorrect focusing	Adjusting the focus system as in the present manual
2. The edge of the image is unfocused.	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
3. Bright patches appear on the image.	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

## **10.0 RECOVERY AND RECYCLING**

Art.13 Dlsg 25 july 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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## **New York Microscope Company Inc**