OPTIKA

M I C R O S C O P E S
I T A L Y

Ver. 4.0.0

 **OPERATION MANUAL**

 **GUIDA UTENTE**

OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0 WARNINGS	page 4
2.0 COMPONENTS	page 6
3.0 INSTALLATION	page 7
4.0 ADJUSTING THE MICROSCOPE	page 11
5.0 PHASE CONTRAST OPERATION	page 15
6.0 MICROPHOTOGRAPHY AND VIDEO	page 16
7.0 TECHNICAL SPECIFICATIONS	page 17
8.0 TROUBLESHOOTING	page 18
9.0 RECOVERY AND RECYCLING	page 21

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

**DANGER**

The instructions on this symbol to avoid possible severe personal injuries.

**WARNING**

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

**WARNING**

Possibility of electric shock.

**HOT!**

Attention: high temperature surfaces. Avoid direct contact.

**NOTE**

Technical notes or usage tips.

SAFETY NOTE

- 1 Do not store the instrument under direct sunlight, high temperature or high humidity, nor in dusty environments. Make sure the stage is plane, horizontal and stable.
- 2 When moving the microscope, please hold the instrument with one hand on the lower side of the eyepiece tube (1), and the other hand on the illumination bracket (2) (Figure 1).
- 3 If water or bacterial cultures are spilled on the stage, objective or viewing tube, pull out the power cord at once, and wipe the microscope clean. Failing to do so may result in damage to the instrument.
- 4 When working, the LED housing on the top of the arm (3) (Figure 1) will become hot. Ensure good ventilation around the LED housing.
- 5 The instrument must be connected to a grounded mains supply.
- 6 Please use the supplied power cord only.

Figure 1

Figure 2


MAINTENANCE

- 1 Use a lint-free cloth to wipe the glass parts. In order to remove fingerprints and oil stains, slightly dampen the cloth with xylene or with a 3:7 mixture of ethanol and ether.
▶ **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- 2 Do not use organic solvents to wipe the non-optical elements. Please use a neutral detergent only.
- 3 During use, if liquids fall on the microscopes, cut off the power at once, and wipe up the moisture.
- 4 Do not disassemble the microscope. Disassembling might degrade the performance of the microscope, and will invalidate the warranty.
- 5 When the objectives are not mounted, please cover the sockets with the dust caps, in order to prevent dust and liquids from entering the microscope.
- 6 When not in use, remember to cover the microscope with the dust hood. Please let the lamp cool down before covering the microscope.

SAFETY SYMBOLS

Symbol	Meaning
	The surface may be hot, do not touch.
	Please read the instruction carefully, before use. Improper operation will result in injuries or malfunction.
	Main switch – ON position.
	Main switch – OFF position.


3.1 INSTALLATION DIAGRAM

The following figure shows the installation sequence of the components.

The number in the figure shows the installation steps.

- ▶ **Before installing, be sure every component is clean. Take care not to scratch parts or surfaces.**
- ▶ **Store the supplied hexagon wrench in a safe place. When changing the components, you will need it again.**

Not included in LT version

3.2 INSTALLATION STEPS

3-2-1 Installing the condenser illumination unit and the LED housing (fig. 3)

1. Insert the condenser illumination unit (1) into the bracket, according to figure .
2. Turn the condenser illumination unit clockwise about 90°, with the “AS” mark of filter holder (3) facing forwards. Align the screw of the condenser illumination unit and the hole of the holder, then screw the bolt in the hole with the supplied hexagon spanner.
3. Insert the connector plug (2) into the connector jack (4).
4. Push the LED housing gently into the holes of the illumination unit.

Figure 3

3-2-2 Installing the objectives fig. 4 & fig. 5

1. Turning the coarse focusing knob (1) till the nosepiece reaches to its lowest position.
 - ▶ **For a safe transport, the nosepiece is placed in the lowest position and the tension adjustment collar (2) is adjusted to the appropriate tension when the microscope leaves the factory.**
2. Screw the lowest magnification objective on to the turret from the right side, then turn the turret clockwise. Mount the other objectives in the same way, following the sequence from low to high.
 - ▶ **Note: the objectives can also be installed through the stage opening.**
 - ▶ Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.
 - ▶ To prevent dust and contamination from entering the microscope, cover all the unused holes with dust caps (3).
 - ▶ When operating, use the low magnification objective (10X) to search and focus the specimen, then switch to higher magnifications.
 - ▶ When switching between objectives, slowly turn the nosepiece until it clicks. The click means that the objective is in the right position, in the center of the light path.


Figure 4


Figure 5

3-2-3 Installing the stage extension and the mechanical stage (fig. 6) (Not included in LT version)

The stage extension can be installed on either side of the stage to enlarge the working surface. The mechanical stage must be installed on the side opposite the extension.

For right-handed operators, the mechanical stage is normally installed on the right side.

1. Installing the stage extension: Screw the bolts (1) on to the extension, then mount the extension from below the stage. Screw it to a firm fit.
2. Installing the mechanical stage: As for the extension, the mechanical stage is fixed by two bolts under the stage.


Figure 6

3-2-4 Installing the stage insert (fig. 7)

1. When using the glass stage (1), make sure that the stage is horizontal.
2. Install the stage insert in the stage opening.

Turn the disk until the V-groove faces the user. This simplifies objective identification.


Figure 7

3-2-5 Installing the eyepieces (fig. 8)

1. Remove the cap of the eyepiece tubes (1).
2. Insert the eyepieces (2) into the tubes.


Figure 8

3-2-6 Installing the color filters (fig. 9)

- ▶ **Be sure the color filter has cooled down completely before you change them. Remove the filter holder (1), then install the color filters (2) you need.**

Mount the color filter flat as shown in (3), verifying that they are not tilted.

- ▶ **If the color filter is tilted or otherwise out of place (4), it may fall.**

The color filters can be stacked in the holder. This allows to install as many filters as needed, as long as the whole thickness is less than 11 mm.

Figure 9

3-2-7 Connecting the power cord (fig.10, 11 & 12)

- ▶ **Do not tension or otherwise stress the power cord. If bent, the wires may easily break.**
 - 1. Turn the main switch (1) to “O”(off) before connecting the power cord.
 - 2. Insert the plug into the power jack (2) of the microscope.
 - 3. Plug the power cord (4) into the mains socket. Check for a safe connection.
- ▶ **Please use the supplied power cord. If lost or damaged, please refer to qualified service.**
 - ▶ **Connect the power cord to a grounded (earthed) power supply only.**


Figure 10

3-2-8 Replacing the fuse (fig. 10 & 11)

Before replacing the fuse, turn the main switch (1) to “O” (off) and unplug the power cord. Rotate the fuse support (3) out of the holder using a straight screwdriver. Insert a new fuse in the support, then rotate the support back into the holder.


Figure 11

- ▶ **Fuse rating: see back of the microscope.**


Figure 12

Aperture diaphragm adjustment lever

Color filter holder

Phase centering bolt

Phase contrast slider


Dioptic adjustment ring

Light path selector lever

Tension adjustment collar

Coarse focus knob

Fine focus knob

4-1 MICROSCOPE BASE

4-1-1 Turning on the LED (fig. 13)

Connect the power, turn on the main switch (1).

4-1-2 Adjusting the brightness (fig. 13)

Turn the brightness adjustment knob (2) to increase and decrease the brightness.


Figure 13

4-1-3 Adjusting the tension (fig. 14)

- ▶ **The coarse focusing knob (1) is preadjusted to a tight tension upon leaving the factory.**

If the nosepiece drops down by itself, or the specimen defocuses while adjusting the fine focus knob (3), the coarse focus knob is too loose. Turning the tension adjustment collar (2) in the indicated direction tightens the coarse focus tension (1). Rotate in the opposite direction to decrease the tension.


Figure 14

4-2 STAGE

4-2-1 Setting the specimen (fig. 15 & 16)

- ▶ **For the best image quality, use flasks, Petri dishes and slides with a 1.2 mm thickness.**

Using $\Phi 35$ mm culture dishes:

You can place a $\Phi 35$ mm culture dish on the stage by using the standard center board (1) of the stage.

Using the mechanical stage:

1. When observing 96-well or 24-well micro-titration plates, please fasten them directly with the stage clip (2).
2. When fastening other plates, please use the adaptor plates supplied with the mechanical stage:
 - Terasaki bracket (3) for Terasaki boards
 - Culture dish bracket (4) for $\Phi 35$ mm culture dishes
 - Object slide bracket (5) for object slides and $\Phi 54$ mm culture dishes
3. Turning the X and Y knobs (6,7), move the specimen to the required position. (Movement Range: 120 (width) \times 78 (length) mm).

Figure 15


Figure 16

(Not included in LT version)

4-2-2 Moving the specimen

Move the specimen to the desired position by freehand or by turning the knobs of the mechanical stage.

- ▶ **When switching objectives, take care not to touch the adaptor plates with the objectives, as their weight may damage the front lens.**

4-3 VIEWING TUBE

4-3-1 Dioptric adjustment (fig. 17)

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring (1) to compensate.
- ▶ **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction.**

Figure 17

4-3-2 Adjusting the interpupillary distance (fig. 18)

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- ▶ **The graduation on the interpupillary distance indicator (3), pointed by the spot “.” (2) on the eyepiece holder, shows the distance between the operator's eyes (fig. 17).**

The range of the interpupillary distance is 48-75mm.

Figure 18

4-3-3 Selecting the light path (fig. 19)

Pull the light path selector lever (1) sideways using your thumb, selecting the light path you need.

Figure 19

	Brightness	Application
IN	100% Used for binocular observation	Binocular observation
OUT	20% Used for binocular observation, and 80% used for video or photography	Binocular observation, television, and micrography or video can be operated simultaneously

4-4 ILLUMINATION UNIT

4-4-1 Using color filters (fig. 20)

Selecting the appropriate color filters according your need.

You can stack a group of color filters in the filter holder, if you ensure that they are level and that the whole thickness is less than 11mm.


Figure 20

Color filter	Use
IF550	Single contrast color filter (green) used for phase contrast microscopy
LBD	Color temperature compensation color filter blue used for bright field observation and microphotography


Figure 21

4-4-2 Using the aperture diaphragm (fig. 21)

When in bright field observation, the aperture diaphragm controls the numerical aperture of the illumination system. When the numerical aperture of the objective and the aperture of the illumination system match, the highest resolution is achieved.

- To identify the aperture diaphragm, the eyepiece can be removed or the centering telescope can be used. The aperture diaphragm looks as in fig. 21.

The aperture can be changed by moving the aperture adjustment lever. ((1)is the image of the aperture diaphragm, (2) is the edge of the objective).

- Generally, when observing a fully chromatic specimen, you need to set the size of the condenser to 70-80% of the aperture of the objective. When observing unstained samples (e.g. bacteria), start from 70% and slowly turn the aperture diaphragm lever clockwise.

5-1 IDENTIFYING THE COMPONENTS

5-1-1 Phase contrast objectives (fig. 22)

The magnifications of the optional phase contrast objectives are: 10X and 20X. These objectives are marked "PH".

Mount the objectives on the turret using the same procedures as the standard objectives (3-2-2).

Figure 22

5-1-2 Phase contrast slider (fig. 23)

Adjustable phase slider.

- The light ring is pre-centered when the microscope leaves the factory. It should therefore need no further adjustment. If a recentering is needed, it can be performed via the two side bolts.
- The 10X/20X light ring (1) must be used with the 10X, 20X phase contrast objectives, the 40X light ring (2) with the 40X and the opening (3) is used for bright field.

Figure 23

5-2 SETUP AND USE

5-2-1 Installing the phase contrast slider (fig. 24)

1. Insert the slider into the illumination system, printed face up.
2. Pull the slider into the desired position, to the click stop.
3. When in phase contrast observation, keep the aperture diaphragm adjustment lever on the "O" (open) position.

Figure 24

5-2-2 Centering the ring (fig. 25)

- ▶ **Usually this operation is not needed. If necessary, please proceed with the following steps:**

1. Place a specimen on the stage and focus it.
2. Take out the eyepiece from the tube without the dioptic adjustment, and replace it with the centering telescope (CT).
3. Check that the phase ring and the objective correspond, and that both are steadily set on a click stop.
4. Use the CT to focus on the light ring's image (5) and the phase contrast ring's image (2). If the light ring's image is not sharp, adjust the CT's eyepiece until you can see a clear image of the light ring (2).
5. Adjust the bolts of the two centering holes in the phase contrast slider using a screwdriver until the light ring center and the phase contrast ring center coincide.
6. The 10X and the 20X phase contrast objectives use the same ring on the phase contrast slider. The coincidence of the light ring center and the phase contrast center must be verified with both objectives.

- ▶ **If the light ring is centered incorrectly, the contrast will be severely impaired.**
- ▶ **The phase ring may need recentering during and after observation of very thick specimens.**
- ▶ **The phase ring may show an apparent misalignment if the cover glass is not flat.**


Figure 25

6-1 VIDEO PORT (VIDEO)

6-1-1 Selecting the light path (fig. 26)

1. To activate the video port, pull out the light path selector lever.
 - ▶ **For observation of dark specimens, first focus, then pull out the lever.**

6-1-2 Installing the video adapter (fig. 27)

1. Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
2. Set up the video adapter on the camera according to its instructions.
3. Install the adapter into the tri-through port, and screw down the bolt (1).

6-1-3 Focus (fig. 27)

During a binocular observation at 20% brightness, look at the image on the video imaging system, refocusing the adapter if necessary.

6-2 MICROPHOTOGRAPHY

6-2-1 Selecting the light path

Set up the light path as in (6-1-1) and in (4-3-3).

6-2-2 Installing the photography adapter (fig. 28)

1. Loosen the locking bolt (1) on the trinocular viewing tube, and take out the dust cap (2).
2. Install the photography adapter into the tri-through port according to its instructions, and screw down the locking bolt (1).
3. Attach the camera ring (if any) to the adapter.
4. Attach the camera to the ring.
 - Warning: for some cameras (mainly reflex) the ring is not included with the microscope, and it should be supplied by the user.
 - For the photography of dark specimens, obscure the eyepieces and the viewfinder with a dark cloth in order to reduce stray light.
 - The camera magnification can be calculated as objective magnification × camera + lens magnification.
 - ▶ **When shooting with a SLR, the mirror movement may cause camera movement. Please lift the mirror, use long exposure times and use an extension cord.**


Figure 26

Figure 27

Figure 28

6-2-3 Focus

During a binocular observation at 20% brightness, look at the image on the video imaging system, refocusing the adapter if necessary.

7.0 TECHNICAL SPECIFICATIONS

7-1 MAIN SPECIFICATIONS

Illumination	Light source type X-LED8 with white 8W LED; light intensity control using a knob on left side of the frame. Color temperature: 6300K LED average life time approx. 50.000h Voltage: 110/240Vac, 50/60Hz, 1A ; Fuse: T500mA 250V Max power required: 13W
Observation Modes	Brightfield, phase contrast
Focusing	Coaxial coarse and fine focusing mechanism (graduated, 0.002mm) with upper stop, to prevent the contact between objective and specimen. Adjustable tension of coarse focusing knob.
Stage	Fixed stage, dimensions 250x230 mm. Mechanical stage mountable on the right side of the stage, X-Y translation range 114x81 mm, with metallic interchangeable inserts for slides, Petri dishes, Terasaki, multi-Well plates, etc. [not included in LT version]. Pair of side extensions to expand the surface of the stage. Glass stage insert with hole for small dimension specimens.
Nosepiece	Quintuple revolving nosepiece, rotation on ball bearings.
Head	Trinocular observation head, inclined 45°. Diopter adjustment on left eyepiece. Interpupillary adjustment 55-75 mm.
Eyepieces	Wide field eyepieces WF10X/22 with field number 22.
Objectives	Infinity corrected optical system IOS (Infinity Optical System). Plan-achromatic LWD objectives infinity corrected, for thickness 1.2 mm, made by following objectives: -) Plan-achromatic IOS LWD 10XPh, N.A. 0.25, W.D. 7.94 mm -) Plan-achromatic IOS LWD 20XPh, N.A. 0.40, W.D. 7.66 mm -) Plan-achromatic IOS LWD 40XPh, N.A. 0.60, W.D. 3.71 mm All objectives are treated with an anti-fungus treatment.
Condenser	LWD condenser, N.A. 0.30, working distance 72 mm. The condenser can be removed to extend the working distance up to 150 mm. Precentered slider with phase rings 10X/20X and 40X
Dimensions	HEIGHT: 485 mm WIDTH: 300 mm DEPTH: 550 mm WEIGHT: 9 kg
Accessories	Hard case [not included in LT version], instruction manual and dust cover included.

8.0 TROUBLESHOOTING

Under certain conditions, some no-fault factors will bring a reversible influence to the instrument's performance. If this happens, please take proper measures according to the follow table. If you can't solve the trouble as indicated, please contact the sales department of our Company.

PROBLEM	REASON	SOLUTION
I. Optical section		
1. The illumination is open, but the field of view is dark.	The plug of the LED holder is not connected to the illumination set	Connect them
	The brightness is too low	Adjust to a proper setting
	Too many colour filters have been stacked	Minimize the number of the filters
2. The edge of the field of view is vignetted or the brightness is asymmetric.	The nosepiece is not in the correct position	Turn the nosepiece to a click stop
	The color filter is partially inserted	Insert the filter to full depth
	The phase contrast slider is not in the proper position	Move the slider to a click stop
3. Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
4. There is an apparent double image.	The size of the aperture diaphragm is too small	Open the aperture diaphragm
5. Poor image quality: The image is not sharp The contrast is not high The details are not clear The phase contrast is low.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The aperture diaphragm in the view of field is opened too much or too little	Adjust the aperture diaphragm
	The lenses (condenser, objective, eyepieces are culture dish) is dirty	Thoroughly clean all the optical system
	In phase contrast observation, the bottom thickness of the sample is more than 1.2mm	Use a sample holder whose bottom thickness is less than 1.2mm
	A bright field objective is used for phase contrast observation	Switch to a phase contrast objective
	The condenser ring is not aligned with the objective phase ring	Adjust the condenser ring to match the objective phase ring
	The light ring and/or the phase contrast ring is not centered	Adjust the bolts to center them
	The objective used is not compatible with the phase ring	Please use a compatible objective
	The phase contrast depends on the sample position	The sample holder is not flat. Move the sample around until a compatible area is found.

6. One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical section		
1. The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
2. The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section		
1. The LED doesn't turn on.	No power supply	Check the power cord connection
2. The brightness is not enough	The brightness adjustment is low	Adjust the brightness
3. The light blinks	The power cord is poorly connected	Check the power cord
IV. Viewing tube assembly		
The field of view of the two eyes is different.	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and video		
1. The image is unfocused.	Incorrect focusing	Adjusting the focus system as in the present manual
2. The edge of the image is unfocused.	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
3. Bright patches appear on the image.	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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New York Microscope Company Inc