

3075-GEM ZOOM STEREO GEM MICROSCOPE

**Designed
For the
Discerning
Gemologist**

ADVANCED
OPTICAL
SYSTEM
FOR EXCELLENT
FLATNESS AND
TREMENDOUS
DEPTH
OF FIELD

EXCEPTIONAL
CLARITY
AND
DETAIL
FOR
VIEWING
DIAMONDS
AND
PRECIOUS
STONES



Delivering quality through vision

3075-GEM ZOOM STEREO GEM MICROSCOPE

Designed for the discerning gemologist, the 3075-GEM incorporates our most popular zoom stereo microscope delivering excellent flatness and tremendous depth of field via our most advanced Greenough optical system ever.

Incorporating the latest ergonomic, optical and mechanical designs, the 3075-GEM Zoom Stereo Microscope delivers unprecedented clarity, detail and accurate color rendition for the examination of diamonds and precious stones.



Accessories shown left:

- ① Inclusion pointer (optional)
- ② Ring holder (optional)
- ③ Wire stone holder (included)

The included wire Stone holder swings out for easier changing of gemstones and rotates for viewing at any angle.

An optional inclusion pointer is also available.



SPECIFICATIONS

Optical System

Featuring a magnification range of 10x—67.5x with the included WF 15x/15mm eyepieces

Illumination

The natural daylight LED illumination system provides dark-field illumination and reflected (top) illumination from above.

The adjustable darkfield condenser with iris allows for precise control of light.

Accessories

Steel wire stone holder allows for multi-angled viewing of various sized gems is included; optional ring holder or inclusion pointer are also available and sold separately

Warranty

Eyeguards, dust cover and manual are also included
5 year limited warranty for parts and labor, and electronic components

Design, features and specifications are subject to change without notice.

The optional ring holder also pivots & rotates for viewing at any angle.



ISO 9001

Design and production adheres to ISO 9001 international quality standard.

ISO 14001

Design and production meets the Requirements of international standard ISO 14001 for environmental management.



New York
Microscope
Company

Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com