

# ACCU-CAM 300 Eyepiece Camera

The ACCU-CAM 300 Eyepiece Camera is the perfect solution for those on a tight budget who need to capture images and video!




CAT# AU-300-EP

- 3.0 MP CMOS eyepiece camera
- Sensor size: 1/2.7"
- USB 2.0 direct connect to your PC!
- Includes our new, feature-rich AccuView software for image analysis, annotation, measurement, white balance & more!
- Captures digital still and video
- Windows compatible
- Easy to use
- Excellent color reproduction
- Fits into a 23.2mm I.D. eye tube or photo tube
- Optional adapter (CAT #5036) available for 30mm I.D. eyepieces
- What's in the box: ACCU-CAM 300 eyepiece camera, USB 2.0 mini cable, AccuView software CD



*Delivering quality through vision*



New York  
**Microscope Company**  
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
Tel: (877) 877-7274 | Fax: (516) 801-2046  
Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)