

# Slide Warmer

The LW Scientific Slide Warmer is designed to be used for a broad range of applications. Our thick metal platform paired with our proprietary heating system ensures stable temperatures and even heating, eliminating the hot and cold spots common to other models. The 12v electronics make it portable when the optional battery cords are purchased – perfect for use by medical and veterinary professionals in the field. Economical and dependable, the LWS Slide Warmer is ideal for use in fertility, cytology, biology, tissue mounting, histology, as well as other applications.

## FEATURES

- Digital temperature display programmable in F° or C°
- PID control ensures accuracy
- Suction cup feet
- Durable metal construction
- Compact, space-saving design
- Thick metal platform
- 14 slide capacity



Designed, Engineered & Assembled Proudly in America

## SPECIFICATIONS

- Temperature range:
  - Ambient +2°F up to 150°F
  - Preset at 98.6°F (37°C)
- Accuracy: +/- 0.3°F

## WEIGHT AND DIMS

- Width: 7.5 in (19 cm)
- Depth: 7.7 in (19.5 cm)
- Height: 4 in (10.7cm)
- Weight: 2.35 lbs (1 kg)

### Platform:

- Width: 7.7" (19.5 cm)
- Depth: 6.9" (17.5 cm)

## ADVANTAGES

- Insulated body
- Compact Design

## Slide Warmer

Model #	Description
SWL-14PL-77DP	Slide Warmer, digital, 14 slide capacity

- Portable – can be operated by 12V battery power
- Even heat distribution
- LED temperature display

## POWER / CONTROLS

- Lighted power switch
- Auto-switching CE/cUL approved power adapter
- Output: 12V8A
- Input: 100-240VAC 50/60Hz

