

# User-friendly, Affordable

## Binocular Stereo Microscope



# EM-20 SERIES

- Robust metal body with smooth rack & pinion  
Fresh design, Ergonomic focus knobs
- Light weight, compact body, easy to carry with hand grip
- Stylish, low-profile light base
- LED features a long life with low power consumption
- Transmitted and incident LED illumination equipped with variable intensity
- Equipped with Rechargeable batteries, cordless  
It enables use of the instrument without an outlet.
- EM-22 has 1X, 2X, 3X Turret
- EM-23 has 1X, 2X, 4X Turret
- Turret Style Objective Changer



MODEL	EM-22	EM-23
Magnification	10X, 20X, 30X	10X, 20X, 40X
Objectives	Turret 1X/2X/3X	Turret 1X/2X/4X
Eyepieces	Paired WF10X/20mm	
Viewing head	Fixed Binocular 45 degree inclined Dioptric adjustment on left eyetube, Interpupillary 55-75mm	
Stand	Rack & pinion stand with Hand grip, Charging function with rechargeable batteries	
Focusing block	Viewing head in up-and-down motion	
Working distance	60 mm	
Illumination	LED incident light/transmitted light, intensity control	
Accessories	Paired eyeshields, AC adapter, dust cover	

USA Branch: URL: <http://www.meijitechno.com>  
 Meiji Techno America  
 5895 Rue Ferrari, San Jose, Ca. 95138  
 Phone: 408-226-3454  
 Fax: 408-226-0900  
 Email: [info@meijitechno.com](mailto:info@meijitechno.com)



Manufacturer: URL: <http://www.meijitechno.com>  
 322-1, Chikumazawa, Miyoshi-machi,  
 Iruma-gun, Saitama 354-0043, Japan  
 Phone : +81-49-259-0111  
 Fax : +81-49-259-0113  
 Email : [meiji@meijitechno.co.jp](mailto:meiji@meijitechno.co.jp)

Óäðã´ ¸áÁ´K

 **New York  
Microscope  
Company**  
 Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
 Tel: (877) 877-7274 | Fax: (516) 801-2046  
 Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)