

# User-friendly, Affordable

## Stereo Microscope

# EM-30 SERIES

- Robust metal body with rack & pinion  
Fresh design, one-point blue parts on white color body
- Light weight, compact body, easy to carry with hand grip
- Stylish, low-profile light base
- Long-life and electricity saving LED light equipped with all models  
Transmitted and incident illumination equipped shown below.



• 2 year Warranty

**EM-32 Model**



### EM-30 Series Ordering Information

MODEL	EM-30	EM-31	EM-32	EM-33
Magnification	10X, 30X	20X, 40X	7X - 45X	
Eyepieces	Paired HWF10X/20mm			
Objectives	Turret 1X/3X	Turret 2X/4X	Zoom 0.7X-4.5X	
Viewing head	Binocular 45 degree inclined			Trino. 45 deg.inclined
Diopter	Dioptric adjustment on both eyetubes, Interpupillary 55-75mm			
Stand	Rack & pinion stand with Hand grip			
Focusing block	Coarse adjustment with tension control			
Working distance	100 mm			
Illumination	3W LED incident light/transmitted light, intensity control,			
Accessories	Paired eyeshields, AC electric cord, dust cover			

**USA Branch:** URL: <http://www.meijitechno.com>  
**Meiji Techno America**  
 5895 Rue Ferrari, San Jose, Ca. 95138  
 Phone: 408-226-3454  
 Fax: 408-226-0900  
 Email: [info@meijitechno.com](mailto:info@meijitechno.com)



**Manufacturer**  
 322-1, Chikumazawa, Miyoshi-machi,  
 Iruma-gun, Saitama 354-0043, Japan  
 Phone : +81-49-259-0111  
 Fax : +81-49-259-0113  
 Email : [meiji@meijitechno.co.jp](mailto:meiji@meijitechno.co.jp)  
 URL: <http://www.meijitechno.com>

Meiji Techno

 **New York  
Microscope  
Company**  
 Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
 Tel: (877) 877-7274 | Fax: (516) 801-2046  
 Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)