

Entry level Student Quality

EDUCATIONAL MICROSCOPE

MT-10 SERIES



MEIJI TECHNO

HIGHLIGHTS

- Ergonomic Hand Grip for portability for ease of storage
- Low positioned Ergonomic Coaxial Focus Control
- Abbe condenser, N.A. 1.25 with iris diaphragm lever for height adjustability
- Revolving Quadruple forward facing nosepiece on ball bearings
- Achromatic DIN Semi Plan Standard 4x, 10X, 40X objective
Optional DIN Semi Plan 60X and 100X objective is available
- Anti-Fungus treated optical parts for mold prevention
- Digital models available with built in 1.3MP 1/2 chip CMOS USB2 Camera
- NeoLED illumination generating even illumination and low electrical consumption
- 5 year Limited Lifetime Warranty

MT-10 Series Ordering Information



MODEL	HEAD	OBJECTIVE	EYEPIECE	STAGE
MT-10	Monocular	4X,10X,40X	DIN WF10X 18mm w/pointer	Plain stage 123mm X 119mm w/ stage clips
MT-11	Monocular	4X,10X,40X,100X	DIN WF10X 18mm w/pointer	Plain stage 123mm X 119mm w/ stage clips
MT-14	Binocular	4X,10X,40X,100X	WF10X 18mm & WF10X w/pointer	Plain stage 123mm x 119mm w/ mechanical stage
MT-16	Digital Monocular	4X,10X,40X,100X	DIN WF10X 18mm w/pointer	Plain stage 123 X 119mm w/ stage clips

USA Branch: URL: <http://www.meijitechno.com>

Meiji Techno America
5895 Rue Ferrari, San Jose, Ca. 95138
Phone: 408-226-3454
Fax: 408-226-0900
Email: info@meijitechno.com



MEIJI TECHNO

Manufacturer: URL: <http://www.meijitechno.com>

322-1, Chikumazawa, Miyoshi-machi,
Iruma-gun, Saitama 354-0043, Japan
Phone : +81-49-259-0111
Fax : +81-49-259-0113
Email : meiji@meijitechno.co.jp

Ōadā eāĀK

**New York
Microscope
Company**
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com