

Advanced Academic observation

EDUCATIONAL MICROSCOPE

MT-30 SERIES



MEIJI TECHNO HIGHLIGHTS

- Ergonomic Hand Grip for portability for ease of storage
- Low positioned Ergonomic Coaxial Focus Control
- Abbe condenser, N.A. 1.25 with iris diaphragm lever for height adjustability
- Revolving Quadruple reverse nosepiece on ball bearings
- Achromatic DIN Semi Plan Standard 4x, 10X, 40X objective
Optional DIN Semi Plan 60X and 100X objective is available
- Anti-Fungus treated optical parts for mold prevention
- Stage lock safety stop to prevent slide glass breakage
- Digital models available with built in 3.2MP 1/3 chip CMOS USB2 Camera
- NeoLED illumination generating even illumination and low electrical consumption
- 2 year Limited Lifetime Warranty



M20 Series Ordering Information

MODEL	HEAD	OBJECTIVE	EYEPIECE	STAGE
MT-30	Binocular	4X,10X,40X,100X	WF10X 18mm & WF10X w/pointer	Plain stage w/mechanical stage 130mm x 130mm
MT-31	Digital Bino.	4X,10X,40X,100X	WF10X 18mm & WF10X w/pointer	Plain stage w/mechanical stage 130mm x 130mm

USA Branch: URL: <http://www.meijitechno.com>

Meiji Techno America

5895 Rue Ferrari, San Jose, Ca. 95138

Phone: 408-226-3454

Fax: 408-226-0900

Emai: info@meijitechno.com



MEIJI TECHNO

Manufacturer: URL: <http://www.meijitechno.com>

322-1, Chikumazawa, Miyoshi-machi,

Iruma-gun, Saitama 354-0043, Japan

Phone : +81-49-259-0111

Fax : +81-49-259-0113

Email : meiji@meijitechno.co.jp

Öä d ä ~ e ä Å ^ K



New York
Microscope
Company
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801

Tel: (877) 877-7274 | Fax: (516) 801-2046

Email: Info@nyscopes.com

www.microscopeinternational.com