

# Moticam **X**

Welcome to the Future of Wireless Digital Microscopy



MICROSCOPES & DIGITAL IMAGING PRODUCTS

For more information, please call 1.800.275.3716



Swift introduces two new unique solutions to turn almost any microscope into a **wireless tablet ready, shareable device**.


## Moticam X



The **Moticam X**, is an “All-In-One-Box” product. What makes the Moticam X a cut above the rest is that it is a high resolution streaming **WiFi camera**. By creating its own wireless network, simply logon with your Android or iOS device and **view, capture, edit live images** from your microscope with our **free App MotiConnect**.

In a case where your school does not allow a download of Apps or you need to use it in a hurry, no problem, as the Moticam X is equipped with its own **on-board software** that **can be controlled simply with a Web Browser** which allows you to capture as well as measure. The Moticam X also works with a WiFi equipped laptop running Windows XP or higher, whereby the camera can be used with standard Motic Images Plus software, enabling you to use this device with your classroom’s Interactive White Board.

## MoticHUB



The **MoticHUB** technology is **integrated** into Motic Images Plus software turning your existing wired Moticam or Digital Microscope (originally shipped with Motic software) into a streaming **WiFi digital solution**.

This **free additional feature** allows you to attach **your digital microscope or camera** to your existing WiFi equipped computer and then stream the image through your existing network for others to pickup on their **WiFi tablets, phones or computers**.



<b>Sensor Type</b>	CMOS
<b>WiFi Resolution</b>	1280x1024 pixels
<b>Chip Size</b>	1/3"
<b>Focusable Lens</b>	12mm
<b>Calibration Slide</b>	Included
<b>Macro Tube</b>	Included
<b>Power Supply</b>	Through provided switching power adapter
<b>Viewing Devices</b>	WiFi enabled PC or Laptop: <i>Motic Images Plus for PC</i> Tablet or Phone Devices: <i>Freely downloadable App MotiConnect</i> Universal Viewer: <i>Almost any HTML 5 supported Web browser</i>

**DISTRIBUTED BY:**

New York Microscope Company Inc.

Phone: (877) 877-7274 | (516) 801-0313

E-mail: [Info@nyscopes.com](mailto:Info@nyscopes.com)

Website: [www.microscopeinternational.com](http://www.microscopeinternational.com)

**MICROSCOPES & DIGITAL IMAGING PRODUCTS**

For more information, please call 1.800.275.3716

