



720p HD C-mount Microscope Camera with Integrated 11,5" HDMI Monitor

Camera specifications	HDMI Camera resolution (MP)	2
	Camera resolution (n° of pixels: W x H)	1280x720
	HDMI Signal output	Yes
	Color / Monochrome	Color
	Sensor Size	1/2.8"
	Sensor Technology	CMOS
	Sensor Type	APTINA CMOS
	Rolling Shutter	Yes
	Image Format	16/9
	Pixel Size (mm)	2.8x2.8
	Frame rate full resolution (fps)	30@1208x720 (HDMI)
	Frame rate other resolution (fps)	30@1920x1080 (Capture)
	Dynamic range (DB)	64
	G sensitivity	510mV at 1/30s
	Dark Signal	0.15mV at 1/30s
	ADC conversion	8 Bit
	Exposure time	0.06 msec - 1.9 sec
	Binning	1x1
	IR filter	Yes
	IR filter range (nm)	380-650 (IR CUT)
	Cooling Temperature	None
	Camera Power	12V 1A
C-Mount	Yes	
Dimensions (mm)	50 x 50 x 65	
Weight (Kg)	0.165	

Camera accessories	HDMI cable (2m)	Yes
	Wireless Mouse	Yes
	16GB SD card	Yes

Environment specifications	Working Temperature	-10°C up to +50°C
	Storage Temperature	-20°C up to +60°C
	Relative Humidity	30% to 80%
	Storage Humidity	10% to 60%

Warranty	Years	2
-----------------	-------	---

HDMI monitor	Size	11,5"
	Power supply	12V / 2,5A
	HDMI cable	50 cm

Óðá dæð ~ e á Á k



100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com