



Model:  
**CL-14**

Typology:  
**ILLUMINATOR**

Description:  
*LED ring light illuminator.*

<b>Cromaticity coordinates (x,y):</b>	0.320, 0.320 (pure white), 6,300 K.
<b>Optical power:</b>	Luminous flux: 240 lumen.  Illuminance (at 10 cm distance): 8,000 lux (on 3 cm diameter).
<b>Power supply:</b>	Brightness adjustment knob. Main power switch 0/1.  External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1,500 mA.
<b>Suitable models:</b>	LAB, SZM and SZN.
<b>Accessories:</b>	Fixing ring for stereo heads.

Optika Microscopes



100 Lauman Lane, Suite A, Hicksville, NY 11801  
Tel: (877) 877-7274 | Fax: (516) 801-2046  
Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)