



Model:
CL-16.1

Typology:
ILLUMINATOR

Description:
LED ring light illuminator.

Cromaticity coordinates (x,y):	0.320, 0.320 (pure white), 6,300 K.
Optical power:	Luminous flux: 500 lumen Illuminance (at 10 cm distance): 6,000 lux (on 10 cm diameter).
Power supply:	Knob for continuous brightness adjustment. Different illumination settings: omnidirectional (ring completely on), oblique (only a quarter or half ring switched on). External power supply provided, with built-in illumination controls. Input: 100/240 Vac 50/60 Hz; output: 12 Vdc 150 mA.
Accessories:	Fixing ring for stereo heads.

Öada~ eãã`K

 **New York
Microscope
Company**
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com