



Model:

CL-30

Typology:

ILLUMINATOR

Description:

Double flexible arm (L = 60 cm) with double LED 2 W.

Chromaticity coordinates (x,y):

0.320, 0.320 (pure white).

Optical power:

Luminous flux: 160 lumen.

Illuminance (at 10 cm distance): 50,000 lux.

Aperture angle:

35°.

Power supply:

External power supply provided.

Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA.

Όραση & Αίτη



100 Lauman Lane, Suite A, Hicksville, NY 11801
 Tel: (877) 877-7274 | Fax: (516) 801-2046
 Email: Info@nyscopes.com
www.microscopeinternational.com