



Model:

CL-41

Typology:

ILLUMINATOR

Description:

Double flexible arm (L = 60 cm) with double X-LED³.

Chromaticity coordinates (x,y): 0.320, 0.320 (pure white).

Electrical power 3.6 W (each arm).

Optical power:
Luminous flux: 400 lumen.
Illuminance (at 10 cm distance): 170,000 lux.

Aperture angle 25°.

Power supply:
Brightness control potentiometer.
External power supply provided.
Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA.

Optika Microscopes



100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com