

**Model:**

Fluorescence attachment with mercury lamp

Description:

Fluorescence with 6-position filterwheel and HBO mercury lamp.

Source	HBO mercury vapour lamp. Power: 100 W. Lifetime: about 300 hours. Knobs for lamp and rear mirror centering.
Filter Set	6 positions fluorescence filter holder: Standard mounted filtersets: - B: EX 460-490, DM 500, EM 510LP; - G: EX 510-550, DM 570, EM 590LP; Optional filtersets: - UV: EX: 330-385 , DM 400 , EM 420LP; - V: EX 400-410, DM 455, EM 455LP. Custom filtersets available on request.
Diaphragms	Aperture and field diaphragms.
Power supply	Included. Voltage: 100/240 Vac, 50/60 Hz Fuse: T3,15 A 250 V. Max power required: 130 W.
Dimensions	620 mm overall length.

Όραση & Ακρόαση



**New York
Microscope
Company**

Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com