

**Model:**

Incident polarizing attachment

Description:*Attachment for incident polarization observation*

Source	High-efficiency X-LED [®] system. Power: 8 W. Lifetime: about 50,000 h.
Filters	Built-in polarizer filter for incident illumination.
Diaphragms	Aperture and field diaphragms.
Power supply	Directly connected to microscope (no external power supply).
Dimensions	600 mm overall length.
Notes	For a complete polarization analysis, use in conjunction with Bertrand Lens attachment.

Όραση και Ακρόαση



Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com