



Model:
Revolving nosepiece

Description:
Revolving nosepiece for RMS objectives with DIC slot

Positions	6
Thread	Standard RMS objectives (20.32 mm thread diameter)
Motorized nosepiece	M-1043 + M-1156 with automatic rotation (in both directions), controlled by buttons, PC software or mouse controller connected to the microscope.
Dimensions	100 mm outer diameter

Όπτικά Οράματα



100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com