



**Model:**  
Revolving nosepiece

**Description:**  
*Revolving nosepiece for polarizing objectives*

<b>Positions</b>	5
<b>Thread</b>	Standard RMS objectives (20.32 mm thread diameter).
<b>Dimensions</b>	85 mm outer diameter.
<b>Notes</b>	Centrable position for each objective (wrench provided).

Öptika Mikroskop



100 Lauman Lane, Suite A, Hicksville, NY 11801  
 Tel: (877) 877-7274 | Fax: (516) 801-2046  
 Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)