



**Model:**  
Revolving nosepiece

**Description:**  
*Revolving nosepiece for incident darkfield objectives*

<b>Positions</b>	5
<b>Thread</b>	For darkfield metallurgical objectives (26 mm thread diameter)
<b>Motorized nosepiece</b>	<b>M-1046 + M-1156</b> with automatic rotation (in both directions), controlled by buttons, PC software or mouse controller connected to the microscope.
<b>Dimensions</b>	100 mm outer diameter
<b>Accessories</b>	3 ring adapters for RMS objectives

Όπτικά Οράκια

 **New York  
Microscope  
Company**  
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
Tel: (877) 877-7274 | Fax: (516) 801-2046  
Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)